

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:** 16/06/2020**Batch No.:****Batch Size:**

|                                       |                                                                                                             |
|---------------------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Generic Name</b>                   | ALLYLESTRENOL TABLETS                                                                                       |
| <b>Label Claim</b>                    | Each Film Coated Tablet contains :<br>Allylestrenol 5 mg<br>Excipients q.s.<br>Colour : Titanium Dioxide IP |
| <b>Mfg. Lic. No.</b>                  | S-MNB/10/93 & S-MB/10/94                                                                                    |
| <b>Pharmacopoeia Status</b>           | IHS                                                                                                         |
| <b>Product Code</b>                   | ALL/001                                                                                                     |
| <b>MFR No.</b>                        | APD/MFR/HT/ALL/001                                                                                          |
| <b>Mfg. Date</b>                      |                                                                                                             |
| <b>Exp. Date</b>                      |                                                                                                             |
| <b>Batch Size</b>                     |                                                                                                             |
| <b>Shelf Life</b>                     | 36 Months                                                                                                   |
| <b>Presentation</b>                   | White, Round, Biconvex, Film Coated Tablets having scored on one side                                       |
| <b>Storage Condition</b>              | Store in a cool, dry, dark place.                                                                           |
| <b>Issued By QA (Sign &amp; Date)</b> |                                                                                                             |
| <b>Issue Date</b>                     |                                                                                                             |
| <b>Date of Commencement</b>           |                                                                                                             |
| <b>Date of Completion</b>             |                                                                                                             |

|                    | Prepared By | Checked By    | Approved By       |
|--------------------|-------------|---------------|-------------------|
| <b>Signature</b>   |             |               |                   |
| <b>Date</b>        | 15/06/2020  | 15/06/2020    | 16/06/2020        |
| <b>Name</b>        | Gulzari Lal | Manjeet Kumar | Namdeo Gath       |
| <b>Designation</b> | Executive   | Sr.Executive  | Jr.Executive      |
| <b>Department</b>  | Production  | Production    | Quality Assurance |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:** 16/08/2020**Batch No.:****Batch Size:****TABLE OF CONTENTS**

| <b>Sr. No.:</b> | <b>Content</b>                                                                             | <b>Page No.</b> |
|-----------------|--------------------------------------------------------------------------------------------|-----------------|
| 1.0             | General Instruction                                                                        | 4               |
| 2.0             | Manufacturing Process Flow Chart                                                           | 5               |
| 3.0             | Master Formula                                                                             | 6               |
| 4.0             | Calculation                                                                                | 7               |
| 5.0             | Dispensing                                                                                 | 8               |
| 5.1             | General Instruction                                                                        | 8               |
| 5.2             | Line Clearance                                                                             | 8-10            |
| 5.3             | Environmental Monitoring                                                                   | 10              |
| 5.4             | RLAF Monitoring                                                                            | 10              |
| 5.5             | Status Labels of Equipment Cleaning                                                        | 11              |
| 5.6             | Dispensing Record                                                                          | 12-14           |
| 6.0             | Manufacturing Record                                                                       | 15              |
| 6.1             | Safety Checklist before Operation of Granulation Area                                      | 15              |
| 6.2             | Details of Equipments /Instruments to be used                                              | 15              |
| 6.3             | Equipments preventive maintenance details                                                  | 16              |
| 6.4             | Equipments/ Instruments cleaning record                                                    | 17              |
| 6.5             | Status labels of equipment's cleaning                                                      | 18              |
| 6.6             | Rinse Water Record (Only in case of Product Change Over)                                   | 19              |
| 6.7             | Line clearance for granulation area                                                        | 20-21           |
| 6.8             | Dispensed Label to be Paste Here                                                           | 22              |
| 6.9             | Dispensed Material Verification                                                            | 23              |
| 6.10            | Sifting                                                                                    | 24              |
| 6.11            | Dry Mixing                                                                                 | 25              |
| 6.12            | Record the weight of sifted and sized granules                                             | 26              |
| 6.13            | Lubricants sifting                                                                         | 27              |
| 6.14            | Blending and Lubrication                                                                   | 27              |
| 6.15            | Lubricated granules weight                                                                 | 28              |
| 6.16            | Blend Sample collection                                                                    | 28              |
| 6.17            | Blend yield Reconciliation                                                                 | 28              |
| 6.18            | Blend release record                                                                       | 29              |
| 7.0             | Compression                                                                                | 29              |
| 7.1             | General Instruction                                                                        | 29              |
| 7.2             | Safety checklist before operation of compression machine                                   | 29              |
| 7.3             | Line Clearance.                                                                            | 30-31           |
| 7.4             | Environmental monitoring of Compression Area.                                              | 32              |
| 7.5             | Compression machine tool Verification.                                                     | 32              |
| 7.6             | Compression Parameters Start Up & In process.                                              | 32              |
| 7.7             | Tablet parameter and in process checks.                                                    | 33              |
| 7.8             | Recording of initial Checks.                                                               | 34              |
| 7.9             | Recording of Initial Checks (at the start up in case of the failure of the initial checks) | 35              |
| 7.10            | Compression Operation Details.                                                             | 36              |
| 7.11            | In-process Checks Quantity & Frequency.                                                    | 36              |
| 7.12            | In process Control Record for Compression.                                                 | 37-40           |
| 7.13            | Core Tablets Weight.                                                                       | 41              |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:**

16/06/2020

**Batch No.:****Batch Size:**

| Sr. No. | Content                                                          | Page No. |
|---------|------------------------------------------------------------------|----------|
| 7.13.3  | Reconciliation of Compressed Tablets.                            | 41       |
| 7.13.4  | Sampling, Analysis and Release of Tablets.                       | 41       |
| 7.13.5  | Transfer above compressed tablets to quarantine                  | 42       |
| 7.13.6  | Issuance of tablets                                              | 42       |
| 8.0     | Coating                                                          | 42       |
| 8.1     | Line Clearance for coating Area                                  | 42-43    |
| 8.2     | Environmental Monitoring of Coating Area                         | 43       |
| 8.3     | Material verification                                            | 44       |
| 8.4     | Sifting                                                          | 44       |
| 8.5     | Solution Preparation                                             | 44-45    |
| 8.6     | Coating Procedure                                                | 45       |
| 8.7     | Time Record for Coating                                          | 46       |
| 8.8     | Coated Tablet Appearance checks                                  | 46       |
| 8.9     | Calculation of % Coating Weight Build up                         | 46       |
| 8.10    | In-Process Specifications                                        | 46       |
| 8.11    | In process Observation                                           | 47       |
| 8.12    | Coated Tablet Weight                                             | 47       |
| 8.13    | Reconciliation of Coated Tablets                                 | 48       |
| 8.14    | Sampling, Analysis and Release of Tablets                        | 48       |
| 8.15    | Transfer above coated tablets to quarantine                      | 48       |
| 8.16    | Issuance of tablets for Inspection                               | 48       |
| 8.17    | Inspection Of Tablets                                            | 49-50    |
| 8.18    | Destruction of Non Recoverable                                   | 50       |
| 9.0     | Accountability                                                   | 51       |
| 10.0    | Signature log                                                    | 51       |
| 11.0    | Incident / deviation received                                    | 51       |
| 12.0    | Break Down Details                                               | 52       |
| 13.0    | Batch manufacturing history sheet                                | 52       |
| 14.0    | Check list for batch manufacturing records & release for packing | 52-53    |
| 15.0    | QC Data review Checklist                                         | 53       |
| 16.0    | Abbreviation                                                     | 53-54    |
| 17.0    | Document revision History                                        | 54       |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS****FOREST ROAD, SOLAN, H.P. (INDIA)****BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:** 16/06/2020**Batch No.:****Batch Size:****1.0 GENERAL INSTRUCTIONS:**

- 1.1 Good manufacturing practices shall be followed during the entire process of manufacturing including sampling and dispensing.
- 1.2 All the containers and equipments used for manufacturing shall be properly cleaned as per the relevant SOP.
- 1.3 All the equipments and containers shall have proper status label as per respective SOP.
- 1.4 All the equipments shall be operated as per the relevant SOP's only.
- 1.5 All the measuring equipments and containers shall be calibrated irrespective of the activity and use within the date of calibration as per the calibration schedule
- 1.6 Weights of all the materials shall be cross check by production personnel against dispensing label before taken up for manufacturing.
- 1.7 Overwriting in BMR and all other documents shall be strictly avoided, and if so, shall be as per APD/QAD
- 1.8 All the operations shall be carried out in clean and orderly manner.
- 1.9 Any deviation in process shall be brought to knowledge of QA and prior approval of QA department shall be taken.
- 1.10 Critical parameters including temperature, humidity and pressure differences shall be checked and monitored APD/QAD.
- 1.11 In process controls shall be carried out throughout manufacturing operations as per relevant instructions defined in this BMR and relevant SOP's.
- 1.12 Ensure that all materials like raw materials, in process materials and finished goods shall be placed in respective areas with proper label.
- 1.13 Attach additional sheets wherever required.
- 1.14 No material shall be placed directly on the floor.



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**


**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**2.0 MANUFACTURING PROCESS FLOW DIAGRAM:**





**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**3.0 MASTER FORMULA:**

| Item Name                           | Specification | Item Code | Mg/Tab  | Overage % | UOM | Std. Qty. for 1,00,000 Tablets |
|-------------------------------------|---------------|-----------|---------|-----------|-----|--------------------------------|
| Allylestrenol                       | IP            | HRMA0001  | 5.000   | 2.0       | Kg  | 0.510                          |
| Lactose Monohydrate                 | IP            | HRME0034  | 115.200 | ----      | Kg  | 11.520                         |
| PVP k-30                            | IP            | HRME0051  | 2.000   | -----     | Kg  | 0.200                          |
| Iso Propyl Alcohol                  | IP            | HRME0081  | 25.000  | -----     | Kg  | 2.500                          |
| Cross carmellose Sodium             | IP            | HRME0014  | 5.000   | -----     | Kg  | 0.500                          |
| Aerosil (Colloidal Silicon Dioxide) | IP            | HRME0002  | 1.000   | -----     | Kg  | 0.100                          |
| Magnesium Stearate                  | IP            | HRME0037  | 1.700   | -----     | Kg  | 0.170                          |
| Protectab HP-1                      | IHS           | HRME0054  | 2.600   | -----     | Kg  | 0.260                          |
| Polysorbate-80 ( Tween-80)          | IP            | HRME0049  | 0.400   | -----     | Kg  | 0.040                          |
| Purified Water                      | IP            | HRME0083  | 23.000  | -----     | Kg  | 2.300                          |



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**



**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date: 16/06/2020**

**Batch No.:**

**Batch Size:**

**4.0 CALCULATION:**

| Ingredient              | Batch Qty. (kg) | Assay                          | API assay on as such basis(A)                                       |
|-------------------------|-----------------|--------------------------------|---------------------------------------------------------------------|
| (a) Allylestrenol       |                 | 100 % assay on as such basis   | $\frac{\text{Assay on dried basis} \times (100 - \text{LOD})}{100}$ |
| (b) Lactose Monohydrate |                 | Final weight to be compensated | =                                                                   |

**The below calculation is to be used when standard quantity of Allylestrenol is to be dispensed is available from one material A.R no.**

Qty. of Allylestrenol is to be dispensed (B) = Std. Qty. X 100 / A = X100/\_\_\_\_\_ = \_\_\_\_\_

Compensated Qty. of Lactose Monohydrate (C) = (a+b) – (B) = \_\_\_\_\_ Kg

**The below table is to be used when standard qty of Allylestrenol is not available from one A.R no.**

| Material A.R. no: _____ | Assay on dried basis (c) | LOD | Assay on as such basis (d) = $\frac{\text{Assay on dried basis} \times (100 - \text{LOD})}{100}$ | Actual quantity available (e) | Qty on 100 % assay basis = (e) X (d)/100 (kg) | Remaining qty to be dispensed (g) = Std. Qty – (f) |
|-------------------------|--------------------------|-----|--------------------------------------------------------------------------------------------------|-------------------------------|-----------------------------------------------|----------------------------------------------------|
|                         |                          |     |                                                                                                  |                               |                                               | (g) = - _____ Kg                                   |
|                         |                          |     |                                                                                                  |                               |                                               | = _____ Kg.                                        |
|                         |                          |     |                                                                                                  |                               |                                               |                                                    |
|                         |                          |     |                                                                                                  | (e) = _____ Kg                | (f) = _____ Kg                                |                                                    |

Next material A.R no: \_\_\_\_\_ Assay on dried basis : \_\_\_\_\_ LOD: \_\_\_\_\_

Assay on as such basis:  $\frac{\text{Assay on dried basis} \times (100 - \text{LOD})}{100} =$  \_\_\_\_\_

Qty. of Allylestrenol is to be dispensed (D) = (g) x 100/ Assay on as such basis = \_\_\_\_\_ Kg

Total Qty. Allylestrenol to be dispensed (H) = (e + D) = \_\_\_\_\_ + \_\_\_\_\_ = \_\_\_\_\_

Compensated Qty. Lactose Monohydrate (I) = (a + b) – H = \_\_\_\_\_ Kg

If assay as such is more than 100% then assay to be consider as 100% for calculation.

\*value depends on LOD and Assay, to be calculated.

**Done By**  
**Production**  
**Sign & Date:**

**Checked By**  
**QA**  
**Sign & Date:**

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:**

16/06/2020

**Batch No.:****Batch Size:****5.0 DISPENSING:****5.1 General Instruction:**

- 5.1 To carry out dispensing of Raw material Follow the respective SOP.
- 5.2 Ensure that the relevant documents shall be updated.
- 5.3 Proper line clearance shall be taken before starting the dispensing activity.
- 5.4 All the activities that are related to Equipment Cleaning, Operations, Material Handling and Process Controls, shall be carried out strictly as per respective Standard Operating Procedure.
- 5.5 Check and verify the Item code No., Material, Quantity & A.R. No. of the Material to be dispensed, is as per Material Requisition Slip prior to moving the material in material entry airlock and performing the dispensing activity.
- 5.6 Verify the Raw Material containers for proper close status before and after performing the dispensing activity.
- 5.7 Take raw materials to dispensing area under RLAF.
- 5.8 Start weighing of active ingredients followed by Excipients.
- 5.9 Weigh and dispense all the materials as per dispensing sheet in double polyethylene bags.
- 5.10 Record the material details in dispensing label and affix the label on the dispensed material and same label shall be placed between two polybag duly approve by Store & QA.
- 5.11 After completion of dispensing activity dispensed material shall transfer in respective production area as per respective SOP.

**5.2 Line Clearance:**To certify the line clearance Follow the SOP No.: **APD/QAD/033.****Date:** \_\_\_\_\_**Time:** \_\_\_\_\_**Previous Product Name:** \_\_\_\_\_**Batch No.:** \_\_\_\_\_

| Sr. No.                                      | Check Points                                                                                         | Observation | Checked By Warehouse Sign & Date | Verified By QA Sign & Date |
|----------------------------------------------|------------------------------------------------------------------------------------------------------|-------------|----------------------------------|----------------------------|
| <b>a] Area General Checks:</b>               |                                                                                                      |             |                                  |                            |
| 1.                                           | Containers/ utensil used for previous batch/product removed from dispensing area                     | Yes/No      |                                  |                            |
| 2.                                           | Containers and vessels in dispensing area shall have cleaned label                                   | Yes/No      |                                  |                            |
| 3.                                           | All documents related to the previous batch/products are removed from dispensing area.               | Yes/No      |                                  |                            |
| 4.                                           | All previous product/batch remnants removed from area                                                | Yes/No      |                                  |                            |
| 5.                                           | Check the Waste bins, it should be clean & empty                                                     | Yes/No      |                                  |                            |
| 6.                                           | Check gowning of the personnel working in the area, gowning shall be done as per respective area SOP | Yes/No      |                                  |                            |
| 7.                                           | Check training records of the personnel working in the area                                          | Yes/No      |                                  |                            |
| <b>b] Area Environmental Control checks:</b> |                                                                                                      |             |                                  |                            |
| 8.                                           | Pressure Difference in Dispensing area<br><b>Limit(NLT 8 Pascal)</b>                                 |             |                                  |                            |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:**

16/06/2020

**Batch No.:****Batch Size:**

| Sr. No.                                     | Check Points                                                   | Observation               | Checked By Warehouse Sign & Date | Verified By QA Sign & Date |
|---------------------------------------------|----------------------------------------------------------------|---------------------------|----------------------------------|----------------------------|
| 9.                                          | Relative Humidity in Dispensing area<br><b>Limit (NMT 55%)</b> |                           |                                  |                            |
| 10.                                         | Temperature of Dispensing area<br><b>Limit (NMT 25°C)</b>      |                           |                                  |                            |
| 11.                                         | Temperature/ humidity log book                                 | Checked/Not checked       |                                  |                            |
| <b>c] Area Cleanliness checks:</b>          |                                                                |                           |                                  |                            |
| 12.                                         | Cleanliness status of floor, corners and walls                 | Cleaned/ not cleaned      |                                  |                            |
| 13.                                         | Area cleanliness record                                        | Checked/not checked       |                                  |                            |
| 14.                                         | Area cleaned by                                                |                           |                                  |                            |
| <b>d] General Checks:</b>                   |                                                                |                           |                                  |                            |
| 15.                                         | Cleanliness status of scoops and other accessories             | Cleaned/ not cleaned      |                                  |                            |
| 16.                                         | Cleanliness status of Balance table.                           | Cleaned/ not cleaned      |                                  |                            |
| <b>e] Equipment Checks:</b>                 |                                                                |                           |                                  |                            |
| <b>Equipment: Balances (APD/WH/BAL/003)</b> |                                                                |                           |                                  |                            |
| 17.                                         | Cleanliness status of Balance                                  | Cleaned/ not cleaned      |                                  |                            |
| 18.                                         | Daily Weight checking log book of balances                     | Checked/not checked       |                                  |                            |
| 19.                                         | Cleanliness beneath the balance                                | Cleaned/ not cleaned      |                                  |                            |
| 20.                                         | Calibration status of Balance                                  | Calibrated/not calibrated |                                  |                            |
| <b>Equipment: Balance (APD/WH/BAL/002)</b>  |                                                                |                           |                                  |                            |
| 21.                                         | Cleanliness status of Balance                                  | Cleaned/ not cleaned      |                                  |                            |
| 22.                                         | Daily Weight checking log book of balances                     | Checked/not checked       |                                  |                            |
| 23.                                         | Cleanliness beneath the balance                                | Cleaned/ not cleaned      |                                  |                            |
| 24.                                         | Calibration status of Balance                                  | Calibrated/not calibrated |                                  |                            |
| <b>Equipment: Balance (APD/WH/BAL/012)</b>  |                                                                |                           |                                  |                            |
| 25.                                         | Cleanliness status of Balance                                  | Cleaned/ not cleaned      |                                  |                            |
| 26.                                         | Daily Weight checking log book of balances                     | Checked/not checked       |                                  |                            |
| 27.                                         | Cleanliness beneath the balance                                | Cleaned/ not              |                                  |                            |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:** 16/06/2020**Batch No.:****Batch Size:**

| Sr. No.                                         | Check Points                                                                       | Observation               | Checked By<br>Warehouse<br>Sign & Date | Verified<br>By QA<br>Sign & Date |
|-------------------------------------------------|------------------------------------------------------------------------------------|---------------------------|----------------------------------------|----------------------------------|
|                                                 |                                                                                    | cleaned                   |                                        |                                  |
| 28.                                             | Calibration status of Balance                                                      | Calibrated/not calibrated |                                        |                                  |
| <b>Equipment: Dispensing Booth (APD/WH/002)</b> |                                                                                    |                           |                                        |                                  |
| 29.                                             | Cleanliness status of dispensing booth                                             | Cleaned/ not cleaned      |                                        |                                  |
| 30.                                             | Cleaning log book of Dispensing record                                             | Checked/not checked       |                                        |                                  |
| 31.                                             | Pressure Difference across HEPA filter in Dispensing Booth<br>(Limit: 8- 20 mm WC) |                           |                                        |                                  |
| 32.                                             | Dispensing booth Cleaned by                                                        |                           |                                        |                                  |

**LINE CLEARANCE CERTIFICATE:**

Observations of Area, Equipment and Environmental checks are found satisfactory / not satisfactory. Line clearance given / not given to proceed with planned process.

**Certified by QA (Sign & Date):** \_\_\_\_\_**Time:** \_\_\_\_\_**5.3 Environmental Monitoring:**

Dispensing of Raw Material shall be carried out under following Environmental Conditions & Record in Environment Monitoring Record at the time of start of dispensing, after every Four Hours:

| Date | Time | Area Name | Temperature<br>NMT25°C | Relative<br>Humidity<br>NMT 55 % | Differential<br>Pressure<br>NLT 08 Pa | Checked By<br>Store<br>Sign & Date | Verified By<br>QA<br>Sign & Date |
|------|------|-----------|------------------------|----------------------------------|---------------------------------------|------------------------------------|----------------------------------|
|      |      |           |                        |                                  |                                       |                                    |                                  |
|      |      |           |                        |                                  |                                       |                                    |                                  |
|      |      |           |                        |                                  |                                       |                                    |                                  |

**Remarks (If any):** \_\_\_\_\_**5.4 RLAF Monitoring:**

| Date | Time | RLAF Reading I.D. No.:APD/WH/002 |                      |                                 |                      | Checked<br>By Store<br>Sign &<br>Date | Verified<br>By QA<br>Sign &<br>Date |
|------|------|----------------------------------|----------------------|---------------------------------|----------------------|---------------------------------------|-------------------------------------|
|      |      | Standard Limit Across            |                      | Observation                     |                      |                                       |                                     |
|      |      | MICROVEE<br>(Pre) Filter         | HEPA(Post)<br>Filter | MICROV<br>EE<br>(Pre)<br>Filter | HEPA(Post)<br>Filter |                                       |                                     |
|      |      | 08-20 mm wc                      | 15-45 mm wc          |                                 |                      |                                       |                                     |
|      |      | 08-20 mm wc                      | 15-45 mm wc          |                                 |                      |                                       |                                     |
|      |      | 08-20 mm wc                      | 15-45mm wc           |                                 |                      |                                       |                                     |
|      |      | 08-20 mm wc                      | 15-45 mm wc          |                                 |                      |                                       |                                     |

**Remarks (If any):** \_\_\_\_\_



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:**

16/06/2020

**Batch No.:**

**Batch Size:**



**5.5 Status Labels of Equipment Cleaning:**

|                                                            |                                               |
|------------------------------------------------------------|-----------------------------------------------|
|                                                            |                                               |
| <b>Label Attached By<br/>Warehouse<br/>Sign &amp; Date</b> | <b>Verified By<br/>QA<br/>Sign &amp; Date</b> |



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**5.6 Dispensing Record:**

**Dispensing of Active Pharmaceutical Ingredients**

**Dispensing Started On:** \_\_\_\_\_

**Dispensing Completed On:** \_\_\_\_\_

| Item Code        | Standard Qty For 1.0 lac Tab. (kg) | Actual Qty For lac Tab. (kg) | UOM | AR No. | Bal. ID No. | Actual weight |          |         | Weigh ed by Store | Ckd By Store Sign & Date | Verified By QA Sign & Date |
|------------------|------------------------------------|------------------------------|-----|--------|-------------|---------------|----------|---------|-------------------|--------------------------|----------------------------|
|                  |                                    |                              |     |        |             | Gross Wt.     | Tare Wt. | Net Wt. |                   |                          |                            |
| Allylestrenol IP |                                    |                              |     |        |             |               |          |         |                   |                          |                            |
| HRMA0001         | 0.510                              |                              |     |        |             |               |          |         |                   |                          |                            |
|                  |                                    |                              |     |        |             |               |          |         |                   |                          |                            |
|                  |                                    |                              |     |        |             |               |          |         |                   |                          |                            |

No. Of Packs/Containers: \_\_\_\_\_

Material Issued By : \_\_\_\_\_

Store Officer  
(Sign & Date)

Material Received By : \_\_\_\_\_

Production  
(Sign & Date)

**Dispensing of Excipients**

**Dispensing Started On:** \_\_\_\_\_

**Dispensing Completed On:** \_\_\_\_\_

| Item Code               | Standard Qty For 1.0 lac Tab. (kg) | Actual Qty For lac Tab. (kg) | UOM | AR No. | Bal. ID No. | Actual weight |          |         | Weighed by Store | Ckd By Store Sign & Date | Verified By QA Sign & Date |
|-------------------------|------------------------------------|------------------------------|-----|--------|-------------|---------------|----------|---------|------------------|--------------------------|----------------------------|
|                         |                                    |                              |     |        |             | Gross Wt.     | Tare Wt. | Net Wt. |                  |                          |                            |
| Lactose Monohydrate IP  |                                    |                              |     |        |             |               |          |         |                  |                          |                            |
| HRME0034                | 11.520                             |                              |     |        |             |               |          |         |                  |                          |                            |
|                         |                                    |                              |     |        |             |               |          |         |                  |                          |                            |
|                         |                                    |                              |     |        |             |               |          |         |                  |                          |                            |
| Lactose Monohydrate IP* |                                    |                              |     |        |             |               |          |         |                  |                          |                            |
| HRME0034                | 11.520                             |                              |     |        |             |               |          |         |                  |                          |                            |
|                         |                                    |                              |     |        |             |               |          |         |                  |                          |                            |
|                         |                                    |                              |     |        |             |               |          |         |                  |                          |                            |

No. Of Packs/Containers: \_\_\_\_\_

Material Issued By : \_\_\_\_\_

Store Officer  
(Sign & Date)

Material Received By : \_\_\_\_\_

Production  
(Sign & Date)



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**Dispensing of Binder**

**Dispensing Started On:** \_\_\_\_\_

**Dispensing Completed On:** \_\_\_\_\_

| Item Code             | Standard Qty For 1.0 lac Tab. (kg) | Actual Qty For lac Tab. (kg) | UOM | AR No. | Bal. ID No. | Actual weight |          |         | Weighed by Store | Ckd By Store Sign & Date | Verified By QA Sign & Date |
|-----------------------|------------------------------------|------------------------------|-----|--------|-------------|---------------|----------|---------|------------------|--------------------------|----------------------------|
|                       |                                    |                              |     |        |             | Gross Wt.     | Tare Wt. | Net Wt. |                  |                          |                            |
| PVP k-30 IP           |                                    |                              |     |        |             |               |          |         |                  |                          |                            |
| HRME0051              | 0.200                              |                              |     |        |             |               |          |         |                  |                          |                            |
|                       |                                    |                              |     |        |             |               |          |         |                  |                          |                            |
|                       |                                    |                              |     |        |             |               |          |         |                  |                          |                            |
| Iso Propyl Alcohol IP |                                    |                              |     |        |             |               |          |         |                  |                          |                            |
| HRME0081              | 2.500                              |                              |     |        |             |               |          |         |                  |                          |                            |
|                       |                                    |                              |     |        |             |               |          |         |                  |                          |                            |
|                       |                                    |                              |     |        |             |               |          |         |                  |                          |                            |

**No. Of Packs/Containers:** \_\_\_\_\_

**Material Issued By :** \_\_\_\_\_  
Store Officer  
(Sign & Date)

**Material Received By :** \_\_\_\_\_  
Production  
(Sign & Date)

**Dispensing of Lubricants**

**Dispensing Started On:** \_\_\_\_\_

**Dispensing Completed On:** \_\_\_\_\_

| Item Code                              | Standard Qty For 1.0 lac Tab(kg) | Actual Qty For lac Tab. (kg) | UOM | AR No. | Bal. ID No. | Actual weight |          |         | Weighed by Store | Checked By Store Sign & Date | Verified By QA Sign & Date |
|----------------------------------------|----------------------------------|------------------------------|-----|--------|-------------|---------------|----------|---------|------------------|------------------------------|----------------------------|
|                                        |                                  |                              |     |        |             | Gross Wt.     | Tare Wt. | Net Wt. |                  |                              |                            |
| Cross carmellose Sodium IP             |                                  |                              |     |        |             |               |          |         |                  |                              |                            |
| HRME0014                               | 0.500                            |                              |     |        |             |               |          |         |                  |                              |                            |
|                                        |                                  |                              |     |        |             |               |          |         |                  |                              |                            |
|                                        |                                  |                              |     |        |             |               |          |         |                  |                              |                            |
| Aerosil (Colloidal Silicon Dioxide) IP |                                  |                              |     |        |             |               |          |         |                  |                              |                            |
| HRME0002                               | 0.100                            |                              |     |        |             |               |          |         |                  |                              |                            |
|                                        |                                  |                              |     |        |             |               |          |         |                  |                              |                            |
|                                        |                                  |                              |     |        |             |               |          |         |                  |                              |                            |
| Magnesium Stearate IP                  |                                  |                              |     |        |             |               |          |         |                  |                              |                            |
| HRME0037                               | 0.170                            |                              |     |        |             |               |          |         |                  |                              |                            |
|                                        |                                  |                              |     |        |             |               |          |         |                  |                              |                            |
|                                        |                                  |                              |     |        |             |               |          |         |                  |                              |                            |

**No. Of Packs/Containers:** \_\_\_\_\_

**Material Issued By :** \_\_\_\_\_  
Store Officer  
(Sign & Date)

**Material Received By :** \_\_\_\_\_  
Production  
(Sign & Date)

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:** 16/06/2020**Batch No.:****Batch Size:****Coating Material**

| Item Code                  | Standard Qty For 1.0 lac Tab (kg) | Actual Qty For lac Tab. (kg) | UOM | AR No. | Bal. ID No. | Actual weight |          |         | Weighed by Warehouse | Ckd By Warehouse Sign & Date | Verified By QA Sign & Date |
|----------------------------|-----------------------------------|------------------------------|-----|--------|-------------|---------------|----------|---------|----------------------|------------------------------|----------------------------|
|                            |                                   |                              |     |        |             | Gross Wt.     | Tare Wt. | Net Wt. |                      |                              |                            |
| Protactab HP- 1            |                                   |                              |     |        |             |               |          |         |                      |                              |                            |
| HRM E0054                  | 0.260                             |                              |     |        |             |               |          |         |                      |                              |                            |
|                            |                                   |                              |     |        |             |               |          |         |                      |                              |                            |
| Polysorbate-80 (Tween -80) |                                   |                              |     |        |             |               |          |         |                      |                              |                            |
| HRM E0049                  | 0.040                             |                              |     |        |             |               |          |         |                      |                              |                            |
|                            |                                   |                              |     |        |             |               |          |         |                      |                              |                            |
| Purified Water             |                                   |                              |     |        |             |               |          |         |                      |                              |                            |
| HRME 0083                  | 2.300                             |                              |     |        |             |               |          |         |                      |                              |                            |
|                            |                                   |                              |     |        |             |               |          |         |                      |                              |                            |

No. Of Packs/Containers: \_\_\_\_\_

Material Issued By : \_\_\_\_\_

Warehouse Officer  
(Sign & Date)

Material Received By : \_\_\_\_\_

Production  
(Sign & Date)



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name:** ALLYLESTRENOL TABLETS

**BMR No.:** APD/BMR/HT/ALL/001

**Revision No.:** 05

**BMR Supersedes No.:** 04

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**6.0 MANUFACTURING RECORD:**

**6.1 Safety Checklist before Operation of Granulation Area:**

**Area:** \_\_\_\_\_


**Date:** \_\_\_\_\_

| Sr. No. | Check Point                                                                        | Observation              | Checked By<br>Production<br>Sign & Date | Verified By<br>QA<br>Sign & Date |
|---------|------------------------------------------------------------------------------------|--------------------------|-----------------------------------------|----------------------------------|
| 1.      | Is any gap between clamp and cover of lighting fixture observed?                   | Complies/Not<br>Complies |                                         |                                  |
| 2.      | Is any electrical cable joint observed?                                            | Complies/Not<br>Complies |                                         |                                  |
| 3.      | Is <b>EMERGENCY STOP</b> switch working (all equipments) can stop instantaneously? | Complies/Not<br>Complies |                                         |                                  |

**Remark (If any):** \_\_\_\_\_

**6.2 Details of Equipments / Instruments to be used:**

| S. No. | Equipment Name                | Equipment Id   | Capacity   | Reference SOP No. |
|--------|-------------------------------|----------------|------------|-------------------|
| 1.     | Balance                       | APD/HT/BAL/001 | 50 Kg      | APD/QAD/033       |
| 2.     | Vibro Sifter                  | APD/HT/003     | 24"        | APD/PRD/002       |
| 3.     | Mass Mixer                    | APD/HT/001     | 50 Kg      | APD/PRD/038       |
| 4.     | Multi mill                    | APD/HT/002     | NA         | APD/PRD/004       |
| 5.     | Tray Dryer                    | APD/HT/004     | 24 Trays   | APD/PRD/037       |
| 6.     | Double Cone Blender           | APD/HT/005     | 50 Kg      | APD/PRD/003       |
| 7.     | Compression Machine           | APD/HT/006     | 30 Station | APD/PGT/006       |
| 8.     | Deduster                      | APD/HT/010     | NA         | APD/PRD/007       |
| 9.     | Colloidal mill                | APD/HT/015     | NA         | APD/PRD/014       |
| 10.    | Coating Pan                   | APD/HT/016     | 36"        | APD/PRD/010       |
| 11.    | Spray Gun                     | APD/HT/021     | 0.5Litre   | APD/PRD/010       |
| 12.    | Hardness Tester               | APD/HT/011     | NA         | APD/QAD/079       |
| 13.    | Disintegration Test Apparatus | APD/HT/013     | NA         | APD/QAD/076       |
| 14.    | Friability Apparatus          | APD/HT/014     | NA         | APD/QAD/074       |
| 15.    | Moisture Analyzer             | APD/HT/030     | NA         | APD/QAD/075       |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:**

16/06/2020

**Batch No.:****Batch Size:****6.3 Equipments Preventive Maintenance Details:**

| S. No. | Equipment Name                | Equipment ID | Detail of Preventive Maintenance |        | Checked By<br>Production<br>Sign & Date | Verified By<br>QA<br>Sign & Date | Remarks<br>(If any) |
|--------|-------------------------------|--------------|----------------------------------|--------|-----------------------------------------|----------------------------------|---------------------|
|        |                               |              | Done on                          | Due on |                                         |                                  |                     |
| 1.     | Vibro Sifter                  | APD/HT/003   |                                  |        |                                         |                                  |                     |
| 2.     | Mass Mixer                    | APD/HT/001   |                                  |        |                                         |                                  |                     |
| 3.     | Multi mill                    | APD/HT/002   |                                  |        |                                         |                                  |                     |
| 4.     | Tray Dryer                    | APD/HT/004   |                                  |        |                                         |                                  |                     |
| 5.     | Double Cone Blender           | APD/HT/005   |                                  |        |                                         |                                  |                     |
| 6.     | Compression Machine           | APD/HT/006   |                                  |        |                                         |                                  |                     |
| 7.     | Deduster                      | APD/HT/010   |                                  |        |                                         |                                  |                     |
| 8.     | Vernier Caliper               | APD/HT/012   |                                  |        |                                         |                                  |                     |
| 9.     | Deduster                      | APD/HT/010   |                                  |        |                                         |                                  |                     |
| 10.    | Colloidal mill                | APD/HT/015   |                                  |        |                                         |                                  |                     |
| 11.    | Coating Pan                   | APD/HT/016   |                                  |        |                                         |                                  |                     |
| 12.    | Spray Gun                     | APD/HT/021   |                                  |        |                                         |                                  |                     |
| 13.    | Disintegration Test Apparatus | APD/HT/013   |                                  |        |                                         |                                  |                     |
| 14.    | Friability Apparatus          | APD/HT/014   |                                  |        |                                         |                                  |                     |
| 15.    | Moisture Analyzer             | APD/HT/030   |                                  |        |                                         |                                  |                     |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:** 16/06/2020**Batch No.:****Batch Size:****6.4 Equipments/Instruments Cleaning Record:**

| S. No. | Machinery/ Equipment      | Capacity      | Equipments ID No. | Previous Product Name | B. No. | Cleaning Status |         | Checked By<br>Producti<br>on Sign &<br>Date | Verifie<br>d By<br>QA Sign &<br>Date |
|--------|---------------------------|---------------|-------------------|-----------------------|--------|-----------------|---------|---------------------------------------------|--------------------------------------|
|        |                           |               |                   |                       |        | Clean On        | Done By |                                             |                                      |
| 1.     | Balance                   | 50 Kg         | APD/HT/BAL/001    |                       |        |                 |         |                                             |                                      |
| 2.     | Vibro Sifter              | 24"           | APD/HT/003        |                       |        |                 |         |                                             |                                      |
| 3.     | Mass Mixer                | 50 Kg         | APD/HT/001        |                       |        |                 |         |                                             |                                      |
| 4.     | Multi mill                | NA            | APD/HT/002        |                       |        |                 |         |                                             |                                      |
| 5.     | Tray Dryer                | 24<br>Trayes  | APD/HT/004        |                       |        |                 |         |                                             |                                      |
| 6.     | Double<br>Cone<br>Blender | 50 Kg         | APD/HT/005        |                       |        |                 |         |                                             |                                      |
| 7.     | Compression<br>Machine    | 30<br>Station | APD/HT/006        |                       |        |                 |         |                                             |                                      |
| 8.     | Deduster                  | NA            | APD/HT/010        |                       |        |                 |         |                                             |                                      |
| 9.     | Vernier<br>Caliper        | NA            | APD/HT/012        |                       |        |                 |         |                                             |                                      |
| 10.    | Colloidal<br>mill         | NA            | APD/HT/015        |                       |        |                 |         |                                             |                                      |
| 11.    | Coating Pan               | 36"           | APD/HT/016        |                       |        |                 |         |                                             |                                      |
| 12.    | Spray Gun                 | 0.5Ltr.       | APD/HT/021        |                       |        |                 |         |                                             |                                      |
| 13.    | Hardness<br>Tester        | NA            | APD/HT/011        |                       |        |                 |         |                                             |                                      |
| 14.    | Disintegratio<br>n Test   | NA            | APD/HT/013        |                       |        |                 |         |                                             |                                      |
| 15.    | Friability<br>Apparatus   | NA            | APD/HT/014        |                       |        |                 |         |                                             |                                      |
| 16.    | Moisture<br>Analyzer      | NA            | APD/HT/030        |                       |        |                 |         |                                             |                                      |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS****FOREST ROAD, SOLAN, H.P. (INDIA)****BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:**

16/06/2020

**Batch No.:****Batch Size:****6.5 Status Labels of Equipment Cleaning:**

**Label Attached By:**  
**Production**  
**Sign & Date**

**Verified By:**  
**QA**  
**Sign & Date**



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date: 16/06/2020**

**Batch No.:**

**Batch Size:**



**6.6 Rinse Water Record (Only in case of Product Change Over):**

| Sr. No. | Equipment Name | Previous Product Name | Batch No. | Cleaning Status (OK/Not OK) | Rinse Water A.R No. | Attached/Not Attached | Checked By             |                |
|---------|----------------|-----------------------|-----------|-----------------------------|---------------------|-----------------------|------------------------|----------------|
|         |                |                       |           |                             |                     |                       | Production Sign & Date | QA Sign & Date |
| 1.      |                |                       |           |                             |                     |                       |                        |                |
| 2.      |                |                       |           |                             |                     |                       |                        |                |
| 3.      |                |                       |           |                             |                     |                       |                        |                |
| 4.      |                |                       |           |                             |                     |                       |                        |                |
| 5.      |                |                       |           |                             |                     |                       |                        |                |
| 6.      |                |                       |           |                             |                     |                       |                        |                |
| 7.      |                |                       |           |                             |                     |                       |                        |                |
| 8.      |                |                       |           |                             |                     |                       |                        |                |
| 9.      |                |                       |           |                             |                     |                       |                        |                |
| 10.     |                |                       |           |                             |                     |                       |                        |                |
| 11.     |                |                       |           |                             |                     |                       |                        |                |

\*Attach the report of Rinse Water Analysis



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**6.7 Line Clearance for Granulation area:**

Before proceeding to the line clearance, ensure that BMR and other related is updated till last stage of processing. To certify the line clearance Follow the SOP No.: APD/QAD/033.

**Date:** \_\_\_\_\_

**Time:** \_\_\_\_\_

**Previous Product Name:** \_\_\_\_\_

**Batch No.:** \_\_\_\_\_

| Sr. No.                                       | Check Points                                                                    | Observation           | Checked By<br>Production<br>Sign & Date | Verified By<br>QA<br>Sign & Date |
|-----------------------------------------------|---------------------------------------------------------------------------------|-----------------------|-----------------------------------------|----------------------------------|
| a]                                            | Area General Checks:                                                            |                       |                                         |                                  |
| 1.                                            | Containers / utensil used for previous batch/product removed from area          | Yes/No                |                                         |                                  |
| 2.                                            | Check the gowning of the working personnel's                                    | Yes/No                |                                         |                                  |
| 3.                                            | Check the training of the working personnel's                                   | Yes/No                |                                         |                                  |
| 4.                                            | Verify Cleaned status label on all equipments for the batch/ product processing | Yes/No                |                                         |                                  |
| 5.                                            | All documents related to the previous batch / product is removed from area.     | Yes/No                |                                         |                                  |
| 6.                                            | All previous product/batch remnants removed from area                           | Yes/No                |                                         |                                  |
| 7.                                            | Check the Waste bins, it should be clean & Empty                                | Yes/No                |                                         |                                  |
| b]                                            | Area Environmental Control checks:                                              |                       |                                         |                                  |
| 8.                                            | Pressure Difference in Granulation area (Limit NLT 08 PASCAL)                   |                       |                                         |                                  |
| 9.                                            | Relative Humidity in Granulation area (Limit NMT 55%)                           |                       |                                         |                                  |
| 10.                                           | Temperature in Granulation area (Limit 25°C)                                    |                       |                                         |                                  |
| 11.                                           | Temperature/Relative humidity log book                                          | Checked / Not checked |                                         |                                  |
| c]                                            | Area Cleanliness checks:                                                        |                       |                                         |                                  |
| 12.                                           | Cleanliness status of floor, corners and walls                                  | Cleaned / Not Cleaned |                                         |                                  |
| 13.                                           | Area cleanliness record                                                         | Cleaned / Not cleaned |                                         |                                  |
| 14.                                           | Area cleaned by                                                                 |                       |                                         |                                  |
| d]                                            | Equipment Checks:                                                               |                       |                                         |                                  |
| Equipment: Vibro sifter (APD/HT/001) & Sieves |                                                                                 |                       |                                         |                                  |
| 15.                                           | Cleanliness status of gasket and clamp                                          | Cleaned/ Not cleaned  |                                         |                                  |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:**

16/06/2020

**Batch No.:****Batch Size:**

| Sr. No.                                                  | Check Points                                             | Observation             | Checked By<br>Production<br>Sign & Date | Verified By<br>QA<br>Sign & Date |
|----------------------------------------------------------|----------------------------------------------------------|-------------------------|-----------------------------------------|----------------------------------|
| 16.                                                      | Cleanliness status of discharge chute                    | Cleaned/<br>Not cleaned |                                         |                                  |
| 17.                                                      | Cleanliness status of Sieve                              | Cleaned/<br>Not cleaned |                                         |                                  |
| 18.                                                      | Integrity of Sieve                                       | Ok/Not Ok               |                                         |                                  |
| <b>Equipment: Mass Mixer (APD/HT/001)</b>                |                                                          |                         |                                         |                                  |
| 19.                                                      | Cleanliness of Shaft and blade                           | Cleaned/<br>Not cleaned |                                         |                                  |
| 20.                                                      | Cleanliness of lid                                       | Cleaned/<br>Not cleaned |                                         |                                  |
| 21.                                                      | Cleanliness of inner bowl surface                        | Cleaned/<br>Not cleaned |                                         |                                  |
| <b>Equipment No: Multi Mill (APD/HT/002)</b>             |                                                          |                         |                                         |                                  |
| 22.                                                      | Blade assembly cleanliness                               | Clean / Not<br>clean    |                                         |                                  |
| 23.                                                      | Integrity of screen                                      | OK/Not OK               |                                         |                                  |
| 24.                                                      | Feed hopper assembly and discharge hopper<br>cleanliness | Clean/ Not<br>clean     |                                         |                                  |
| <b>Equipment: Double Cone Blender 50 KG (APD/HT/004)</b> |                                                          |                         |                                         |                                  |
| 25.                                                      | Cleanliness status of Pillar Blender                     | Cleaned/<br>not cleaned |                                         |                                  |

**LINE CLEARANCE CERTIFICATE:**

Observations of Area, Equipment and Environmental checks are found satisfactory/not satisfactory. Line clearance given / not given to proceed with planned process.

**Certified by QA (Sign & Date):** \_\_\_\_\_**Time:** \_\_\_\_\_



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**6.8 Dispensed Label to be Paste Here:**

**Label Attached By:**  
**Production**  
**Sign & Date**

**Verified By:**  
**QA**  
**Sign & Date**

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:**

16/06/2020

**Batch No.:****Batch Size:****6.9 Dispensed Material Verification**

| Sr. No     | Item code | Material                            | Specification | AR No. | Gross Weight (kg) | Total Quantity (kg) | Checked By Production Sign & Date | Verified By QA Sign & Date |
|------------|-----------|-------------------------------------|---------------|--------|-------------------|---------------------|-----------------------------------|----------------------------|
| DRY MIXING |           |                                     |               |        |                   |                     |                                   |                            |
| 1.         | HRMA0001  | Allylestrenol                       | IP            |        |                   |                     |                                   |                            |
| 2.         | HRME0034  | Lactose Monohydrate                 | IP            |        |                   |                     |                                   |                            |
| 3.         | HRME0051  | PVP k-30                            | IP            |        |                   |                     |                                   |                            |
| 4.         | HRME0081  | Iso Propyl Alcohol                  | IP            |        |                   |                     |                                   |                            |
| 5.         | HRME0014  | Croscarmellose Sodium               | IP            |        |                   |                     |                                   |                            |
| 6.         | HRME0002  | Aerosil (Colloidal Silicon Dioxide) | IP            |        |                   |                     |                                   |                            |
| 7.         | HRME0037  | Magnesium Stearate                  | IP            |        |                   |                     |                                   |                            |

Material Verified By: \_\_\_\_\_  
Production  
(Sign & Date)

Material Verified By: \_\_\_\_\_  
QA  
(Sign & Date)



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**6.10 Sifting:**

**6.10.1** Masks and gloves should be used during sifting.

**6.10.2** Check for the line clearance, Cleanliness status of the area and equipment.

**6.10.3** Check and record the Temperature, Humidity and pressure difference.

**Frequency: At the start of Operation and every four hours thereafter and end of the process.**

| Date | Time | Area Name | Temperature<br>(NMT25°C) | Relative<br>Humidity<br>NMT 55% | Differentia<br>l Pressure<br>NLT 8 Pa | Checked By<br>Production<br>Sign & Date | Verified By<br>QA<br>Sign & Date |
|------|------|-----------|--------------------------|---------------------------------|---------------------------------------|-----------------------------------------|----------------------------------|
|      |      |           |                          |                                 |                                       |                                         |                                  |
|      |      |           |                          |                                 |                                       |                                         |                                  |

**Remarks (If Any):**

**6.10.4** Transfer the dispensed material to manufacturing area from day store area for verification.

**6.10.5** Check the sieve integrity before and after sifting as per SOP No.: APD/PRD/006

**6.10.6** Operate the Sifter as per the SOP APD/PRD/ 002 for sifting.

**Date of Operation:**

**Start Time:**

**End Time:**

| Sr.<br>No. | Material Name          | Specific<br>ation | Total Container<br>wise Qty. Gross<br>wt.<br>(in kg) | Total<br>Net<br>Qty.<br>In Kg | Sieve<br>Size | Sieve Integrity |       | Sifting Time |    | Done By<br>Production<br>Sign &<br>Date | Checked By<br>QA Sign &<br>Date |
|------------|------------------------|-------------------|------------------------------------------------------|-------------------------------|---------------|-----------------|-------|--------------|----|-----------------------------------------|---------------------------------|
|            |                        |                   |                                                      |                               |               | Before          | After | From         | To |                                         |                                 |
| 1.         | Allylestrenol          | IP                |                                                      |                               | 60#           |                 |       |              |    |                                         |                                 |
| 2.         | Lactose<br>Monohydrate | IP                |                                                      |                               | 40 #          |                 |       |              |    |                                         |                                 |

**6.11 Binder Preparation**

**Date of Operation:** **Start Time:** **End Time:**

**6.11.1** Masks and gloves should be used during operation.

**6.11.2** Check for the line clearance, cleanliness of the area and equipment.

**6.11.3** Operate the paste preparation kettle as per the SOP APD/PGT/026

**6.11.4** Take Iso Propyl Alcohol in a Paste Kettle.

**6.11.5** Add polyvinyl pyrrolidone K-30 (PVPK-30) with continuous stirring.

| S. No. | Item Name                               | Quantity | Added By | Checked By Production<br>Sign & Date |
|--------|-----------------------------------------|----------|----------|--------------------------------------|
| 1.     | Polyvinyl pyrrolidone K-30<br>(PVPK-30) |          |          |                                      |
| 2.     | Iso Propyl Alcohol                      |          |          |                                      |

**Appearance of Binder:-**

**Checked By:**



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**



**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**

**6.12 Dry Mixing:**

**Date of Operation:** \_\_\_\_\_

**6.12.1** Masks and gloves should be used during operation.

**6.12.2** Check for the line clearance, cleanliness of the area and equipment.

**6.12.3** Operate the Mass Mixer as per the SOP APD/PRD/038 for dry mixing.

| S. No. | Item Name           | Quantity kg | Added into Mass Mixer By | Checked By Production Sign & Date |
|--------|---------------------|-------------|--------------------------|-----------------------------------|
| 1.     | Allylestrenol       |             |                          |                                   |
| 2.     | Lactose Monohydrate |             |                          |                                   |

**6.12.4** Mix the materials in Mass Mixer for 30 min.

| Date | Time Start | Time End | Done By | Checked By Prod. Sign & Date |
|------|------------|----------|---------|------------------------------|
|      |            |          |         |                              |

**6.13 Wet Granulation:**

**Date of Operation:** \_\_\_\_\_

**6.13.1** Masks and hand gloves should be used during operation.

**6.13.2** Check for the line clearance, cleanliness of the area and equipment.

**6.13.3** Check and record the Temperature and Humidity as per SOP APD/QAD/037.

**6.13.4** Operate the Mass mixer as per the SOP APD/PRD/038 for wet granulation as per following parameter.

**6.13.5** Add the binder solution to Mass mixer and mix it for 10 mins.

**6.13.6** After achieve the required granule stop the Mass Mixer.

**6.13.7** Parameters for wet granulation:

| Parameter                        | STD Time          | Start Time | End Time | Total Time in minute | Done By | Checked By Prod. Sign |
|----------------------------------|-------------------|------------|----------|----------------------|---------|-----------------------|
| Addition of binder solution      | Within 10 Minutes |            |          |                      |         |                       |
| Granulation                      | Within 12 Minutes |            |          |                      |         |                       |
| Unload wet mass into FBD trolley | 2-3 Minutes       |            |          |                      |         |                       |

**6.14 Drying:**

**Date of Operation:** \_\_\_\_\_

**6.14.1** Masks and gloves should be used during operation.

**6.14.2** Check for the line clearance, cleanliness of the area and equipment.

**6.14.3** Fix the Trays with wet granules, Switch on the Heater and switch on the drier and dry the granules for 45-60 minutes at an inlet temperature of 50-60° C with an outlet temperature 40±5°c.

**6.14.4** Put off the drier, remove the trays of tray dryer and rake the materials.

| Date | Status | Time |     | Duration in mins | Temp (° C) |        | Done By | Checked By Sign & Date |
|------|--------|------|-----|------------------|------------|--------|---------|------------------------|
|      |        | On   | Off |                  | Inlet      | Outlet |         |                        |
|      | Drying |      |     |                  |            |        |         |                        |



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**6.15 Sizing of Semi-dried granule :**

**Date of Operation:** \_\_\_\_\_

**6.15.1** Masks and gloves should be used during operation.

**6.15.2** Milled over size granules through 2.5 mm Multi mill screen

| Date | Screen Size | Screen No. | Screen Integrity |           | Start Time | End Time | Total Time Min. | Done By | Checked By Sign & Date |
|------|-------------|------------|------------------|-----------|------------|----------|-----------------|---------|------------------------|
|      |             |            | Before Use       | After Use |            |          |                 |         |                        |
|      | 2.5mm       |            |                  |           |            |          |                 |         |                        |

**6.16 Final Drying:**

**Date of Operation:** \_\_\_\_\_

**6.16.1** Drying the milled granules for 10 minutes at inlet temperature 50-60°C and an outlet temperature 40°C. check the LOD of the granules. If required dry until loss on drying is between 0.2 – 1.8 %w/w.

**6.16.2** After regular interval, turn the mass of trays upside down with the help of SS scoop.

| Date | Status       | Time |     | Duration in mins | Temp (°C) |        | Done By | Checked By Sign & Date | LOD % | Done By | Checked By Sign & Date |
|------|--------------|------|-----|------------------|-----------|--------|---------|------------------------|-------|---------|------------------------|
|      |              | On   | Off |                  | Inlet     | Outlet |         |                        |       |         |                        |
|      | Drying Final |      |     |                  |           |        |         |                        |       |         |                        |

**6.17 Final Sizing & Sifting :-**

**6.17.1** Sift final dried granules through sifter using Mesh No.24#.

| Date | Sieve Size | Sieve ID No. | Screen Integrity |           | Start Time | End Time | Total Time Min. | Done By | Checked By Sign & Date |
|------|------------|--------------|------------------|-----------|------------|----------|-----------------|---------|------------------------|
|      |            |              | Before Use       | After Use |            |          |                 |         |                        |
|      | 24 #       |              |                  |           |            |          |                 |         |                        |

**6.17.2** Mill the remaining residue through 1.5mm Multimill screen.

| Date | Screen Size | Screen No. | Screen Integrity |           | Start Time | End Time | Total Time Min. | Done By | Checked By Sign & Date |
|------|-------------|------------|------------------|-----------|------------|----------|-----------------|---------|------------------------|
|      |             |            | Before Use       | After Use |            |          |                 |         |                        |
|      | 1.5 mm      |            |                  |           |            |          |                 |         |                        |

**6.11 Record the weight of sifted Material :-**

| S. No.                | Gross Weight (In Kg) | Tare Weight (In Kg) | Net Weight (In Kg) | Done By Sign & Date | Checked By Sign & Date |
|-----------------------|----------------------|---------------------|--------------------|---------------------|------------------------|
| 1                     |                      |                     |                    |                     |                        |
| 2                     |                      |                     |                    |                     |                        |
| 3                     |                      |                     |                    |                     |                        |
| TOTAL NET WEIGHT (KG) |                      |                     |                    |                     |                        |

**Theoretical Weight of Sifted Material: -** \_\_\_\_\_ **kg**

**Remarks (IfAny):** \_\_\_\_\_



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**6.12 Lubricants sifting.**

**Operation Start Time:**

**Operation End Time:**

| Sr. No | Material Name             | Specification | Total Bag/ Container wise Qty. Gross wt. (in kg) | Total Net Qty. In Kg | Sieve Size | Sieve Integrity |       | Sifting Time |    | Done By Prod. Sign & Date | Checked By QA Sign & Date |
|--------|---------------------------|---------------|--------------------------------------------------|----------------------|------------|-----------------|-------|--------------|----|---------------------------|---------------------------|
|        |                           |               |                                                  |                      |            | Before          | After | From         | To |                           |                           |
| 1.     | Cross Carmellose Sodium   | IP            |                                                  |                      | 40#        |                 |       |              |    |                           |                           |
| 2.     | Magnesium Stearate        | IP            |                                                  |                      | 40#        |                 |       |              |    |                           |                           |
| 3.     | Colloidal silicon dioxide | IP            |                                                  |                      | 40#        |                 |       |              |    |                           |                           |

**6.13 Blending/Lubrication:**

**6.14.1** Operate the Blender as per the SOP No. **APD/PRD/003**.

**6.14.2** Masks and gloves should be used during operation.

**6.14.3** Check for the line clearance, cleanliness of the area and equipment.

**6.14.4** Charge the above mixture of Norethisterone & Lactose granules in blender.

**6.14.5** Charge the lubrication and lock the blender properly.

**6.14.6** Mix the blend in blender for 40 min with following parameter.

| Date | RPM of Double cone blender |          | Time Start | Time End | Done By | Checked By Production Date & Sign |
|------|----------------------------|----------|------------|----------|---------|-----------------------------------|
|      | Std                        | Observed |            |          |         |                                   |
|      | 10                         |          |            |          |         |                                   |

**6.14.7** Now Stop the blender & Charge the sifted Magnesium Stearate, Mix the blend in blender for 05 min with following parameter.

| Date | RPM of Double cone blender |          | Time Start | Time End | Done By | Checked By Production Date & Sign |
|------|----------------------------|----------|------------|----------|---------|-----------------------------------|
|      | Std                        | Observed |            |          |         |                                   |
|      | 10                         |          |            |          |         |                                   |

**6.14.8** After Completion of blending activity transfer the lubricated blend in Double Poly Bag lined HDPE drum/ bin, affix Status and Under Test label as per SOP No.: **APD/QAD/011** and transfer to Granules quarantine area, weigh and record the details.



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**6.14 Lubricated Granules Weight:**

| Drum No.               | Clean Status<br>OK/Not ok | Gross Wt.<br>(Kg) | Tare Wt.<br>(Kg) | Net Wt.<br>(Kg) | Done By<br>Production | Checked By<br>Production<br>Sign & Date | Remarks<br>(If any) |
|------------------------|---------------------------|-------------------|------------------|-----------------|-----------------------|-----------------------------------------|---------------------|
| 1.                     |                           |                   |                  |                 |                       |                                         |                     |
| 2.                     |                           |                   |                  |                 |                       |                                         |                     |
| 3.                     |                           |                   |                  |                 |                       |                                         |                     |
| 4.                     |                           |                   |                  |                 |                       |                                         |                     |
| Total No. of Container |                           | Total weight      |                  |                 |                       |                                         |                     |

**Theoretical weight of Blend:** \_\_\_\_\_ kg.

**Limit: To be ascertained**

**Remarks (If Any):** \_\_\_\_\_

**Physical Parameters for lubricated Granules:-**

| Parameters           | Limits | Unit | Observation | Instrument Code | Checked by (Sign & Date) |
|----------------------|--------|------|-------------|-----------------|--------------------------|
| L.O.D. of Lubricated |        | %w/w |             |                 |                          |

**6.15 Blend Sample Collection Record:**

Send the sample intimation to QA for sampling.

After sampling, QA shall send the sample along-with intimation to QC for analysis.

Sampling of Lubricated granules shall be performed as per APD/QAD/

| Intimation No. | Intimated By Production<br>Date | Sampled By QA<br>Date | Quantity Sampled |
|----------------|---------------------------------|-----------------------|------------------|
|                |                                 |                       |                  |
|                |                                 |                       |                  |

**6.16 Yield Calculation: (Lubricated Blend)**

|    |                                            |                       |         |                               |
|----|--------------------------------------------|-----------------------|---------|-------------------------------|
| 1. | Theoretical weight of blend                |                       |         | kg.                           |
| 2. | Actual Weight of blend                     |                       |         | Kg                            |
| 3. | In-process samples                         |                       |         |                               |
|    | A                                          | IPQA Samples (if any) | .....kg |                               |
|    | B                                          | QC Samples            | .....kg |                               |
|    | C                                          | Validation Samples    | .....kg |                               |
| 4. | Total samples: (A + B + C)                 |                       |         | kg                            |
| 5. | Total: (2 + 4)                             |                       |         | Kg                            |
| 6. | Non-Recoverable                            |                       |         |                               |
| 7. | Yield $\{(5/1) \times 100\} =$ _____ X 100 |                       |         | _____%                        |
|    | (Record the yield)                         |                       |         | <b>Limit: To be ascertain</b> |

**Remarks (If Any):** \_\_\_\_\_

**Done By Production**  
(Sign & Date)

**Checked By Production**  
(Sign & Date)

**Verified By QA**  
(Sign & Date)

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:**

16/06/2020

**Batch No.:****Batch Size:****6.17 Blend Release Record:**

(To be filled by QA Officer /Executive)

Date of Release: \_\_\_\_\_

Sample Intimation No.: \_\_\_\_\_

AR No.: \_\_\_\_\_

Result: **Complies / does not Complies****Test Report attached/Not Attached****Checked By QA****(Sign & Date)** \_\_\_\_\_**7.0 COMPRESSION:****7.1 General Instruction:****7.1.1** Gowning should be as per respective SOP.**7.1.2** Check for the line clearance, Cleanliness status of the area and equipment.**7.1.3** After receiving release of lubricated granules from Quality Assurance, proceed for the compression.**7.1.4** Carry out the Initial Checks of Compression Machine setting and physical parameters of tablets and record.**7.1.5** Mask and gloves should be used throughout compression.**7.1.6** Counter check weights of total finished granules.**7.1.7** Set the punch & dies tools & operate the compression machine, Deduster as per SOP  
**APD/PRD/007****7.1.8** In case if compression activity is carried out on next day i.e. process is started after shutdown; additional pages for compression activity (Including Safety Checklist, Line clearance, Environmental Monitoring of Compression Area, Compression Machine Tool Verification, Compression Parameters Start Up & In-process, Compression Operation Details.**7.2 Safety Checklist (Before Operation of Compression Machine):****Area:** \_\_\_\_\_ **Equipment ID No:** \_\_\_\_\_ **Date:** \_\_\_\_\_

| Sr. No. | Check Point                                                                             | Observation              | Checked By<br>Production<br>Date & Sign | Verified By<br>QA<br>Date & Sign |
|---------|-----------------------------------------------------------------------------------------|--------------------------|-----------------------------------------|----------------------------------|
| 1.      | Does the Door interlocks (limit switches) working properly?                             | Complies/Not<br>Complies |                                         |                                  |
| 2.      | Are all required Personnel Protective Equipments available?                             | Complies/Not<br>Complies |                                         |                                  |
| 3.      | Are Electrical cables in good condition and free from joints?                           | Complies/Not<br>Complies |                                         |                                  |
| 4.      | Is <b>EMERGENCY STOP</b> switch working (Compression machine) can stop instantaneously? | Complies/Not<br>Complies |                                         |                                  |

**Remarks (If Any):**



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**7.3 Line Clearance:**

To certify the line clearance Follow the SOP No.: APD/QAD/034.


**Date:** \_\_\_\_\_

**Time:** \_\_\_\_\_

**Previous Product Name:** \_\_\_\_\_

**Batch No.:** \_\_\_\_\_

| Sr. No                                             | Check Points                                                                    | Observation             | Checked By |    |
|----------------------------------------------------|---------------------------------------------------------------------------------|-------------------------|------------|----|
|                                                    |                                                                                 |                         | Production | QA |
| a) Area General Checks:                            |                                                                                 |                         |            |    |
| 1.                                                 | Containers/ Utensil materials used for previous batch/product removed from area | Yes/No                  |            |    |
| 2.                                                 | Cleaned status label on all equipments for the batch/product                    | Yes/No                  |            |    |
| 3.                                                 | All documents related to the previous batch / products are removed from area.   | Yes/No                  |            |    |
| 4.                                                 | BMR till last processed stage is updated                                        | Yes/No                  |            |    |
| 5.                                                 | Check the Waste bins; it should be clean & empty.                               | Yes/No                  |            |    |
| b) Area Environmental Control Checks:              |                                                                                 |                         |            |    |
| 6.                                                 | Pressure Difference in Compression area (Limit: NLT8Pascal)                     |                         |            |    |
| 7.                                                 | Relative Humidity in Compression area (Limit: NMT 55%)                          |                         |            |    |
| 8.                                                 | Temperature of Compression area (Limit: NMT 25 °C)                              |                         |            |    |
| 9.                                                 | Temperature/Relative humidity log book                                          | Checked/<br>Not checked |            |    |
| c) Area Cleanliness checks:                        |                                                                                 |                         |            |    |
| 10.                                                | Cleanliness status of floor, corners and walls                                  | Cleaned/ Not<br>cleaned |            |    |
| 11.                                                | Area cleanliness record                                                         | Cleaned/Not<br>cleaned  |            |    |
| 12.                                                | Area cleaned by                                                                 |                         |            |    |
| d) Equipment No.: Compression Machine (APD/HT/006) |                                                                                 |                         |            |    |
| 13.                                                | Cleanliness status of Feeder assembly                                           | Cleaned/Not<br>cleaned  |            |    |
| 14.                                                | Cleanliness status of Hopper assembly                                           | Cleaned/Not<br>cleaned  |            |    |
| 15.                                                | Cleanliness status of Turret, hole and scrapper                                 | Cleaned/Not<br>cleaned  |            |    |
| 16.                                                | Cleanliness status of roller                                                    | Cleaned/Not<br>cleaned  |            |    |
| 17.                                                | Cleanliness and No. of the Lower cam track                                      | Ok/Not ok               |            |    |
| e) Equipment No.: Deduster: (APD/HT/010)           |                                                                                 |                         |            |    |
| 18.                                                | Tablet Path cleanliness                                                         | Cleaned/Not             |            |    |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:**

16/06/2020

**Batch No.:****Batch Size:**

| Sr. No.                                               | Check Points                                                 | Observation               | Checked By |    |
|-------------------------------------------------------|--------------------------------------------------------------|---------------------------|------------|----|
|                                                       |                                                              |                           | Production | QA |
|                                                       |                                                              | cleaned                   |            |    |
| 19.                                                   | Cleanliness status of Dust collector pipe and assembly       | Cleaned/Not cleaned       |            |    |
| 20.                                                   | Cleanliness of electrical panel                              | Cleaned/Not cleaned       |            |    |
| <b>f] Punch Tool Checks:</b>                          |                                                              |                           |            |    |
| 21.                                                   | Punch tool No. and Identity as per BMR                       | Ok/ Not Ok                |            |    |
| 22.                                                   | Cleaning log book check for defined Punch tool               | Checked/Not checked       |            |    |
| 23.                                                   | Previous Product: Punch tool No. used                        |                           |            |    |
| 24.                                                   | Cleanliness status of Punches                                | Cleaned/Not cleaned       |            |    |
| 25.                                                   | Setting of punches on Compression machine                    | Ok/Not Ok                 |            |    |
| <b>g] Equipment: Balance (BAL:- (APD/HT/BAL/002))</b> |                                                              |                           |            |    |
| 26.                                                   | Cleanliness status of Balance                                | Cleaned/ not cleaned      |            |    |
| 27.                                                   | Daily Weight checking log book of balances                   | Checked/not checked       |            |    |
| 28.                                                   | Calibration status of Balance                                | Calibrated/not calibrated |            |    |
| <b>h] Miscellaneous Equipments:</b>                   |                                                              |                           |            |    |
| 29.                                                   | Calibration status of Hardness Tester (APD/HT/011)           | Calibrated/not calibrated |            |    |
| 30.                                                   | Calibration status of Vernier Caliper (APD/HT/012)           | Calibrated/not calibrated |            |    |
| 31.                                                   | Calibration status of Moisture analyzer balance (APD/HT/030) | Calibrated/not calibrated |            |    |
| 32.                                                   | Calibration status of Disintegration Apparatus (APD/HT/013)  | Calibrated/not calibrated |            |    |
| 33.                                                   | Calibration status of Friability Apparatus (APD/HT/014)      | Calibrated/not calibrated |            |    |

**LINE CLEARANCE CERTIFICATE:**

Observations of Area, Equipment and Environmental checks are found satisfactory/ not satisfactory. Line clearance given/ not given to proceed with planned process.

**Certified by QA (Sign & Date):** \_\_\_\_\_

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:**

16/06/2020

**Batch No.:****Batch Size:****7.4 Environmental Monitoring of Compression Area:**

Check and record the Temperature, Humidity and pressure difference as per SOP No.: APD/QAD/007.

**Frequency:** At the start of operation and thereafter every four hours.

| Date | Time | Area Name | Temp.<br>NMT<br>25°C | Relative<br>Humidity<br>NMT<br>55% | Differential<br>Pressure<br>NLT8Pascal | Checked By<br>Production<br>Date & Sign | Verified By<br>QA<br>Date & Sign |
|------|------|-----------|----------------------|------------------------------------|----------------------------------------|-----------------------------------------|----------------------------------|
|      |      |           |                      |                                    |                                        |                                         |                                  |
|      |      |           |                      |                                    |                                        |                                         |                                  |
|      |      |           |                      |                                    |                                        |                                         |                                  |
|      |      |           |                      |                                    |                                        |                                         |                                  |

**Remarks (If Any):** \_\_\_\_\_**7.5 Compression Machine Tool Verification:****Date:** \_\_\_\_\_**Time:** \_\_\_\_\_

| Particulars              | Dies | Upper Punch               | Lower Punch               | Checked By<br>Production<br>Sign & Date | Verified By<br>QA<br>Sign & Date | Remark<br>(If any) |
|--------------------------|------|---------------------------|---------------------------|-----------------------------------------|----------------------------------|--------------------|
| Punch size               | NA   | 7.10 mm                   | 7.10 mm                   |                                         |                                  |                    |
| Shape<br>Specification   | NA   | White, Round.<br>Biconvex | White, Round.<br>Biconvex |                                         |                                  |                    |
| Marking<br>Specification | NA   | Scored on one<br>Side     | NA                        |                                         |                                  |                    |
| Quantity of<br>Dies      |      |                           |                           |                                         |                                  |                    |
| Quantity of<br>Punches   |      |                           |                           |                                         |                                  |                    |

**7.6 Compression Parameters Start Up & In-process:****7.6.1** Carry out in-process as per SOP No.: APD/QAD/037**7.6.2** Rotate hand wheel clockwise for two rotations and weigh the tablets.**7.6.3** Set the weight of tablet at 130.00 mg by turning weight dossier anticlockwise for increasing weight and clockwise for decreasing weight. Discard the tablets.**7.6.4** Operate machine in auto mode. Check the total weight of 20 tablets.**7.6.5** Set weight by turning the Weight dozer anticlockwise for increasing the weight.**7.6.6** Set weight by turning the Weight dozer clockwise for decreasing the Weight.**7.6.7** Check the weight of 20 tablets & calculate the average weight of tablet.**7.6.8** Check the Individual Weight Variation, Diameter, Thickness And Hardness.**7.6.9** Set thickness and hardness as per limits given below by turning punch penetration handle clockwise for increasing hardness and reducing thickness and by turning punch penetration handle anticlockwise for decreasing hardness and increasing thickness for each side.**7.6.10** Check friability, disintegration time and individual weight variation of tablet. Note down the data as start up checks.



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**


**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date: 16/06/2020**

**Batch No.:**

**Batch Size:**



**7.6.11 Compression Machine Setting Parameters: Equipment ID No.: APD/HT/006**

| M/c Parameters  | Acceptance Criteria                    |
|-----------------|----------------------------------------|
| Compression M/c | 30 Station Single Rotary, (BB tooling) |
| Upper Punch     | 7.1 mm Round, Biconvex                 |
| Lower Punch     | 7.1 mm Round, Biconvex                 |
| Machine Speed   |                                        |

**7.7 Tablet Parameters and in process checks:**

| Tablet Parameters                  | Acceptance Criteria                                       |
|------------------------------------|-----------------------------------------------------------|
| Average wt of tablet               | 130 mg $\pm$ 7.5%                                         |
| Limits of Individual tablet Weight | 120.25 mg - 139.75 mg                                     |
| Thickness                          | 3.2 mm $\pm$ 0.2 mm (3.0mm-3.4mm)                         |
| Limits of Group Tablets Weight     | 2.600gm $\pm$ 2% (2.548gm -2.652 gm)                      |
| Disintegration Time                | N.M.T. 15 Min.                                            |
| Hardness (Kg/ cm <sup>2</sup> )    | NLT 3.0 Kg/cm <sup>2</sup>                                |
| Friability (%)                     | N.M.T. 1.0 % w/w                                          |
| Appearance                         | White, Round, Biconvex, Scored one Side, Uncoated Tablets |

**Formula for Calculating Friability:**

$$\text{Friability} = \frac{(\text{Initial weight} - \text{Weight after rotation})}{\text{Initial weight}} \times 100$$



**Batch Size:**



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD,SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.:05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**7.9 Recording of Initial Checks (at the start up in case of the failure of the initial checks):**

| Parameters                                                                                                   | Limits                                                    | Observations |  |     |      |  |     |
|--------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|--------------|--|-----|------|--|-----|
|                                                                                                              |                                                           | Production   |  |     | IPQA |  |     |
|                                                                                                              | Date & Time                                               |              |  |     |      |  |     |
|                                                                                                              | NA                                                        | LHS          |  | RHS | LHS  |  | RHS |
| Machine Speed                                                                                                | From.....to .....                                         |              |  |     |      |  |     |
| Physical Appearance                                                                                          | White, Round, Biconvex, Scored one Side, Uncoated Tablets |              |  |     |      |  |     |
| Deduster                                                                                                     | OK / Not OK                                               |              |  |     |      |  |     |
| Group Weight of 20 units                                                                                     | 2.600gm $\pm$ 2%<br>(2.548–2.652 gm)                      |              |  |     |      |  |     |
| Individual weight of 30 units                                                                                | 130 mg $\pm$ 7.5%<br>(120.25mg – 139.75 mg)               |              |  |     |      |  |     |
|                                                                                                              |                                                           |              |  |     |      |  |     |
|                                                                                                              |                                                           |              |  |     |      |  |     |
|                                                                                                              |                                                           |              |  |     |      |  |     |
|                                                                                                              |                                                           |              |  |     |      |  |     |
|                                                                                                              |                                                           |              |  |     |      |  |     |
|                                                                                                              |                                                           |              |  |     |      |  |     |
|                                                                                                              |                                                           |              |  |     |      |  |     |
|                                                                                                              |                                                           |              |  |     |      |  |     |
|                                                                                                              |                                                           |              |  |     |      |  |     |
| Thickness of 10 units                                                                                        | 3.2 mm $\pm$ 0.2 mm<br>(3.0mm-3.4mm)                      |              |  |     |      |  |     |
| Hardness of 10 units                                                                                         | NLT 3.0 Kg/cm <sup>2</sup>                                |              |  |     |      |  |     |
| Friability = $\frac{\text{Initial Wt.} - \text{Final Wt.}}{\text{Initial Wt.}} \times 100$<br>(NMT 1.0% w/w) |                                                           |              |  |     |      |  |     |
| Disintegration Time                                                                                          | NMT 15 min                                                |              |  |     |      |  |     |
| Officer/Executive Sign & Date                                                                                |                                                           |              |  |     |      |  |     |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:** 16/06/2020**Batch No.:****Batch Size:****7.10 Compression Operation Details:**

| Date        | Compression Started at | Compression Completed at | Operated By | Checked By<br>Production<br>Sign & Date |
|-------------|------------------------|--------------------------|-------------|-----------------------------------------|
|             |                        |                          |             |                                         |
|             |                        |                          |             |                                         |
|             |                        |                          |             |                                         |
| Total Time: |                        |                          |             |                                         |

**7.11 In-process Checks Quantity & Frequency:**

| Parameters                              | Frequency                                                           |
|-----------------------------------------|---------------------------------------------------------------------|
| Appearance of units of one round        | <b>Production:</b> Every One hour<br><br><b>QA:</b> Every Two hours |
| Individual Weight Variation of 30 units |                                                                     |
| Group Wt. of 20 units                   |                                                                     |
| Thickness of 10 units                   |                                                                     |
| Hardness of 10 units                    |                                                                     |
| Friability: 50 units                    | <b>Production:</b> Every two hour<br><b>QA:</b> Every Four hours    |
| Disintegration Time of 6 units          |                                                                     |

**Quantity to be taken for in-process check: \_\_\_\_\_ Units (Each Side)**



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**7.12 In process Control Record for Compression:**

| Parameters                                                                                                   | Limits                                                    | Observations |     |     |     |  |  |
|--------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|--------------|-----|-----|-----|--|--|
|                                                                                                              | Date                                                      |              |     |     |     |  |  |
|                                                                                                              | Start Time                                                |              |     |     |     |  |  |
|                                                                                                              | NA                                                        | LHS          | RHS | LHS | RHS |  |  |
| Machine Speed                                                                                                | From.....to .....                                         |              |     |     |     |  |  |
| Physical Appearance                                                                                          | White, Round, Biconvex, Scored one Side, Uncoated Tablets |              |     |     |     |  |  |
| Deduster                                                                                                     | OK / Not OK                                               |              |     |     |     |  |  |
| Group Weight of 20 units                                                                                     | 2.600gm $\pm$ 2%<br>(2.548-2.652 gm)                      |              |     |     |     |  |  |
| Individual weight of 30 units                                                                                | 130 mg $\pm$ 7.5%<br>(120.25mg- 139.75 mg)                |              |     |     |     |  |  |
|                                                                                                              |                                                           |              |     |     |     |  |  |
|                                                                                                              |                                                           |              |     |     |     |  |  |
|                                                                                                              |                                                           |              |     |     |     |  |  |
|                                                                                                              |                                                           |              |     |     |     |  |  |
|                                                                                                              |                                                           |              |     |     |     |  |  |
|                                                                                                              |                                                           |              |     |     |     |  |  |
|                                                                                                              |                                                           |              |     |     |     |  |  |
|                                                                                                              |                                                           |              |     |     |     |  |  |
|                                                                                                              |                                                           |              |     |     |     |  |  |
| Thickness of 10 units                                                                                        | 3.2mm $\pm$ 0.2 mm<br>(3.0mm-3.4mm)                       |              |     |     |     |  |  |
| Hardness of 10 units                                                                                         | NLT 3.0 Kg/cm <sup>2</sup>                                |              |     |     |     |  |  |
| Friability = $\frac{\text{Initial Wt.} - \text{Final Wt.}}{\text{Initial Wt.}} \times 100$<br>(NMT 1.0% w/w) |                                                           |              |     |     |     |  |  |
| Disintegration Time                                                                                          | NMT 15 min                                                |              |     |     |     |  |  |
| Checked By<br>Production/QA<br>Sign and Date                                                                 |                                                           |              |     |     |     |  |  |

**Remark (If any):** \_\_\_\_\_



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**7.12 In process Control Record for Compression:**

| Parameters                                                                                                   | Limits                                                    | Observations |     |     |     |     |     |
|--------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|--------------|-----|-----|-----|-----|-----|
|                                                                                                              | Date                                                      |              |     |     |     |     |     |
|                                                                                                              | Start Time                                                |              |     |     |     |     |     |
|                                                                                                              | NA                                                        | LHS          | RHS | LHS | RHS | LHS | RHS |
| Machine Speed                                                                                                | From.....to .....                                         |              |     |     |     |     |     |
| Physical Appearance                                                                                          | White, Round, Biconvex, Scored one Side, Uncoated Tablets |              |     |     |     |     |     |
| Deduster                                                                                                     | OK / Not OK                                               |              |     |     |     |     |     |
| Group Weight of 20 units                                                                                     | 2.600gm $\pm$ 2%<br>(2.548–2.652 gm)                      |              |     |     |     |     |     |
| Individual weight of 30 units                                                                                | 130 mg $\pm$ 7.5%<br>(120.25mg – 139.75 mg)               |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
| Thickness of 10 units                                                                                        | 3.2mm $\pm$ 0.2 mm<br>(3.0mm–3.4mm)                       |              |     |     |     |     |     |
| Hardness of 10 units                                                                                         | NMT 3.0 kg/cm <sup>2</sup>                                |              |     |     |     |     |     |
| Friability = $\frac{\text{Initial Wt.} - \text{Final Wt.}}{\text{Initial Wt.}} \times 100$<br>(NMT 1.0% w/w) |                                                           |              |     |     |     |     |     |
| Disintegration Time                                                                                          | NMT 15 min                                                |              |     |     |     |     |     |
| Checked By Prod./QA Sign and Date                                                                            |                                                           |              |     |     |     |     |     |

**Remark (If any):** \_\_\_\_\_

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:** 16/06/2020**Batch No.:****Batch Size:****7.12 In process Control Record for Compression:**

| Parameters                                                                                                   | Limits                                                    | Observations |     |     |     |     |     |
|--------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|--------------|-----|-----|-----|-----|-----|
|                                                                                                              | Date                                                      |              |     |     |     |     |     |
|                                                                                                              | Start Time                                                |              |     |     |     |     |     |
|                                                                                                              | NA                                                        | LHS          | RHS | LHS | RHS | LHS | RHS |
| Machine Speed                                                                                                | From.....to .....                                         |              |     |     |     |     |     |
| Physical Appearance                                                                                          | White, Round, Biconvex, Scored one Side, Uncoated Tablets |              |     |     |     |     |     |
| Deduster                                                                                                     | OK / Not OK                                               |              |     |     |     |     |     |
| Group Weight of 20 units                                                                                     | 2.600gm $\pm$ 2%<br>(2.548–2.652 gm)                      |              |     |     |     |     |     |
| Individual weight of 30 units                                                                                | 130 mg $\pm$ 7.5%<br>(120.25mg – 139.75 mg)               |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
| Thickness of 10 units                                                                                        | 3.2mm $\pm$ 0.2 mm<br>(3.0mm-3.4mm)                       |              |     |     |     |     |     |
| Hardness of 10 units                                                                                         | NLT 3.0 Kg/cm <sup>2</sup>                                |              |     |     |     |     |     |
| Friability = $\frac{\text{Initial Wt.} - \text{Final Wt.}}{\text{Initial Wt.}} \times 100$<br>(NMT 1.0% w/w) |                                                           |              |     |     |     |     |     |
| Disintegration Time                                                                                          | NMT 15 min                                                |              |     |     |     |     |     |
| Checked By Prod./QA Sign and Date                                                                            |                                                           |              |     |     |     |     |     |

Remark (If any):

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:** 16/06/2020**Batch No.:****Batch Size:****7.12 In process Control Record for Compression:**

| Parameters                                                                                                   | Limits                                                    | Observations |     |     |     |     |     |
|--------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|--------------|-----|-----|-----|-----|-----|
|                                                                                                              | Date                                                      |              |     |     |     |     |     |
|                                                                                                              | Start Time                                                |              |     |     |     |     |     |
|                                                                                                              | NA                                                        | LHS          | RHS | LHS | RHS | LHS | RHS |
| Machine Speed                                                                                                | From.....to .....                                         |              |     |     |     |     |     |
| Physical Appearance                                                                                          | White, Round, Biconvex, Scored one Side, Uncoated Tablets |              |     |     |     |     |     |
| Deduster                                                                                                     | OK / Not OK                                               |              |     |     |     |     |     |
| Group Weight of 20 units                                                                                     | 2.600gm $\pm$ 2%<br>(2.548–2.652 gm)                      |              |     |     |     |     |     |
| Individual weight of 30 units                                                                                | 130 mg $\pm$ 7.5%<br>(120.25mg – 139.75 mg)               |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
|                                                                                                              |                                                           |              |     |     |     |     |     |
| Thickness of 10 units                                                                                        | 3.2mm $\pm$ 0.2 mm<br>(3.0mm-3.4mm)                       |              |     |     |     |     |     |
| Hardness of 10 units                                                                                         | NLT 3.0 Kg/cm <sup>2</sup>                                |              |     |     |     |     |     |
| Friability = $\frac{\text{Initial Wt.} - \text{Final Wt.}}{\text{Initial Wt.}} \times 100$<br>(NMT 1.0% w/w) |                                                           |              |     |     |     |     |     |
| Disintegration Time                                                                                          | NMT 15 min                                                |              |     |     |     |     |     |
| Checked By Prod./QA Sign and Date                                                                            |                                                           |              |     |     |     |     |     |

**Remark (If any):** \_\_\_\_\_



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001** **Revision No.: 05**

**BMR Supersedes No.: 04** **Effective Date: 16/06/2020**

**Batch No.:** **Batch Size:**



**7.13 Core Tablets Weight:**

**7.13.1 Verify the Balance Calibration Status prior to weighing:**

| Balance ID | Capacity | Calibration Status<br>(Ok/ Not Ok) | Checked By QA<br>Sign & Date | Verified By QA<br>Sign & Date |
|------------|----------|------------------------------------|------------------------------|-------------------------------|
|            |          |                                    |                              |                               |

**7.13.2 Weighing Record:**

| Container No.          | Clean status<br>Ok/Not Ok | Gross Wt.<br>(Kg) | Tare Wt.<br>(Kg) | Net Wt.<br>(Kg) | Remarks<br>(If any) |
|------------------------|---------------------------|-------------------|------------------|-----------------|---------------------|
| 1.                     |                           |                   |                  |                 |                     |
| 2.                     |                           |                   |                  |                 |                     |
| 3.                     |                           |                   |                  |                 |                     |
| 4.                     |                           |                   |                  |                 |                     |
| Total No. of Container |                           |                   | Total weight     |                 |                     |

**Theoretical weight of compressed Tablet:** \_\_\_\_\_ **kg.**

**Limit: To be ascertained**

**Done By Production:** \_\_\_\_\_  
**(Sign & Date)**

**Checked By Q.A.** \_\_\_\_\_  
**(Sign & Date)**

**7.13.3 Reconciliation of Compressed Tablets:**

| Sr. No. | Calculation Points                    | Observed Values | Checked By<br>Production<br>Sign & Date | Verified By<br>QA<br>Sign & Date |
|---------|---------------------------------------|-----------------|-----------------------------------------|----------------------------------|
| 1.      | Theoretical Batch size (Kg) (A)       |                 |                                         |                                  |
| 2.      | Weight of Compressed Tablets (Kg) (B) |                 |                                         |                                  |
| 3.      | In process Sample (a)                 |                 |                                         |                                  |
| 4.      | QC Sample (b)                         |                 |                                         |                                  |
| 5.      | Validation Sample (c)                 |                 |                                         |                                  |
| 6.      | Total Sample {(a)+ (b) + (c) = (C)}   |                 |                                         |                                  |
| 7.      | Actual Yield [(B+C/A) x100]           |                 |                                         |                                  |
| 8.      | Non Recoverable (Kg)                  |                 |                                         |                                  |

**7.13.4 Sampling, Analysis and Release of Tablets**

Send the intimation slip to QA Dept. to draw the sample of Tablets for analysis.

**Sampling, Analysis Details.**

| Intimated By Production<br>Sign & Date | Intimation Received By<br>QA Sign & Date | Sampled By QA<br>Sign, Date & Time | Quantity Sampled |
|----------------------------------------|------------------------------------------|------------------------------------|------------------|
|                                        |                                          |                                    |                  |

The compressed tablets are **Released/ Not Released** for Inspection/ Packing.

**QA Sign & Date:** \_\_\_\_\_

**Time:** \_\_\_\_\_

**Attach the Report from QC:** \_\_\_\_\_

**A.R. No:** \_\_\_\_\_



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date: 16/06/2020**

**Batch No.:**

**Batch Size:**



**7.13.5 Transfer compressed tablets to quarantine.**

| Date of Transfer | Qty. Transferred (in kg) | Transferred By Sign & Date | Verified By QA Sign & Date |
|------------------|--------------------------|----------------------------|----------------------------|
|                  |                          |                            |                            |

**7.13.6 Issuance of tablets:**

| Date of Transfer | Date of Receiving | Qty. Received in kg | Received By Sign & Date | Verified By QA Sign & Date |
|------------------|-------------------|---------------------|-------------------------|----------------------------|
|                  |                   |                     |                         |                            |

**8.0 Coating :**

**8.1 Line Clearance for coating Area:**

To certify the line clearance Follow the SOP No.: **APD/QAD/034**

**Date:** \_\_\_\_\_

**Time:** \_\_\_\_\_

**Previous Product Name:** \_\_\_\_\_

**Batch No.:** \_\_\_\_\_

| S. No. | Check Points | Observation | Checked By             |                |
|--------|--------------|-------------|------------------------|----------------|
|        |              |             | Production Sign & Date | QA Sign & Date |

**a) Area General Checks:**

|    |                                                                                 |        |  |  |
|----|---------------------------------------------------------------------------------|--------|--|--|
| 1. | Containers/ Utensil materials used for previous batch/product removed from area | Yes/No |  |  |
| 2. | Cleaned status label on all equipments for the batch/product                    | Yes/No |  |  |
| 3. | All documents related to the previous batch /Product are removed from area.     | Yes/No |  |  |
| 4. | BMR till last processed stage is updated                                        | Yes/No |  |  |
| 5. | All previous product/batch remnants removed from area                           | Yes/No |  |  |
| 6. | Check the Waste bins, it should be clean & empty                                | Yes/No |  |  |
| 7. | Check gowning of the personnel working in the area                              | Yes/No |  |  |
| 8. | Check training record of the personnel working in the area                      | Yes/No |  |  |

**b) Area Environmental Control checks:**

|     |                                                               |  |  |  |
|-----|---------------------------------------------------------------|--|--|--|
| 9.  | Pressure Difference in Dispensing area (Limit: NLT 08 Pascal) |  |  |  |
| 10. | Relative Humidity in Dispensing area (Limit: NMT 55%)         |  |  |  |



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



| S. No.                                       | Check Points                                             | Observation             | Checked By                |                   |
|----------------------------------------------|----------------------------------------------------------|-------------------------|---------------------------|-------------------|
|                                              |                                                          |                         | Production<br>Sign & Date | QA<br>Sign & Date |
| 11.                                          | Temperature of Dispensing area<br>(Limit: NMT 25°C)      |                         |                           |                   |
| 12.                                          | Temperature/Relative humidity log book                   | Checked/<br>Not checked |                           |                   |
| c) Area Cleanliness checks:                  |                                                          |                         |                           |                   |
| 13.                                          | Cleanliness status of floor, corners and walls           | Checked/<br>Not checked |                           |                   |
| 14.                                          | Area cleanliness record                                  | Checked/<br>Not checked |                           |                   |
| 15.                                          | Area cleaned by                                          |                         |                           |                   |
| d) Equipment No: Manual Coating (APD/GT/009) |                                                          |                         |                           |                   |
| 16.                                          | Check spray gun, dosing delivery line for<br>cleanliness | Checked/<br>Not checked |                           |                   |
| 17.                                          | Solution tank and stirrer cleanliness                    | Checked/<br>Not checked |                           |                   |
| 18.                                          | Check pan and below baffles for the cleanliness          | Checked/Not<br>checked  |                           |                   |

**LINE CLEARANCE CERTIFICATE:**

Observations of Area, Equipment and Environmental checks are found satisfactory/not satisfactory.

Line clearance given/not given to proceed with planned process.

**Certified by QA (Sign & Date):** \_\_\_\_\_

**Time:** \_\_\_\_\_

**8.2 Environmental Monitoring of Coating Area:**

Check and record the Temperature, Humidity and pressure difference as per SOP No.: APD/QAD/

**Frequency:** At the start of operation and thereafter every four hours.

**Area Name:**

| Date | Time | Temperature<br>NMT 25°C | Relative Humidity<br>NMT 55% | Differential Pressure<br>NLT 08 Pascal | Checked By<br>Production<br>Date & Sign | Verified By<br>QA<br>Date & Sign |
|------|------|-------------------------|------------------------------|----------------------------------------|-----------------------------------------|----------------------------------|
|      |      |                         |                              |                                        |                                         |                                  |
|      |      |                         |                              |                                        |                                         |                                  |
|      |      |                         |                              |                                        |                                         |                                  |
|      |      |                         |                              |                                        |                                         |                                  |
|      |      |                         |                              |                                        |                                         |                                  |

**Remarks (If Any):** \_\_\_\_\_



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name:** ALLYLESTRENOL TABLETS

**BMR No.:** APD/BMR/HT/ALL/001

**Revision No.:** 05

**BMR Supersedes No.:** 04

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**8.3 Material verification :**

Transfer the dispensed material to manufacturing area from day store area for verification.

**Start Time:** \_\_\_\_\_

**End Time:** \_\_\_\_\_

| S. No. | Material Name  | Specification | Total Bag/<br>Container wise Qty.<br>Gross wt.<br>(in kg) | Total<br>Net Qty.<br>In Kg | Checked By<br>Production Sign& Date | Verified By<br>QA Sign&<br>Date |
|--------|----------------|---------------|-----------------------------------------------------------|----------------------------|-------------------------------------|---------------------------------|
| 1.     | Protectab HP-1 | IHS           |                                                           |                            |                                     |                                 |
|        |                |               |                                                           |                            |                                     |                                 |
|        |                |               |                                                           |                            |                                     |                                 |
|        |                |               |                                                           |                            |                                     |                                 |
| 2.     | Tween-80       | IP            |                                                           |                            |                                     |                                 |
|        |                |               |                                                           |                            |                                     |                                 |
|        |                |               |                                                           |                            |                                     |                                 |
|        |                |               |                                                           |                            |                                     |                                 |
| 3.     | Purified Water | IHS           |                                                           |                            |                                     |                                 |
|        |                |               |                                                           |                            |                                     |                                 |
|        |                |               |                                                           |                            |                                     |                                 |
|        |                |               |                                                           |                            |                                     |                                 |

**8.4 Sifting :-**

**8.4.1** Transfer the dispensed material to manufacturing area from day store area for verification.

**8.4.2** Check the sieve integrity before and after sifting as per SOP No.: APD/PRD/002

**8.4.3** Operate the Sifter as per the SOP APD/PRD/002 for sifting.

**Date of Operation:** \_\_\_\_\_

**Start Time:** \_\_\_\_\_

**End Time:** \_\_\_\_\_

**8.5 Solution Preparation :-**

**Date of Operation:** \_\_\_\_\_

**Start Time:** \_\_\_\_\_

**End Time:** \_\_\_\_\_


**8.5.1** Take \_\_\_\_\_ kg Purified Water in S.S. Vessel/Container.

**8.5.2** Add Sifted Protectab HP-1 & Tween-80 into above SS Vessel/Container with constant stirring.  
Mix well for **15 min.**

**8.5.3** Rinse it with of Purified Water.



**Batch Size:**



| Date | Time Start | Time End | Done By | Checked By Prod.<br>Sign & Date |
|------|------------|----------|---------|---------------------------------|
|      |            |          |         |                                 |

**8.6.9** Check the Parameter at the time of start and then follow the given frequency.

Page No.: 45 of 54



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**8.7 Time Record for Coating:**

**Date:**

| Loading |    | Jogging |    | Warming |    | Spraying |    | Cooling & Drying |    | Unloading of Tablets |    | Total Coating Time | Operator | Checked By Production Sign & Date |
|---------|----|---------|----|---------|----|----------|----|------------------|----|----------------------|----|--------------------|----------|-----------------------------------|
| From    | To | From    | To | From    | To | From     | To | From             | To | From                 | To |                    |          |                                   |
|         |    |         |    |         |    |          |    |                  |    |                      |    |                    |          |                                   |
|         |    |         |    |         |    |          |    |                  |    |                      |    |                    |          |                                   |

**8.8 Coated Tablet Appearance checks:**

| Appearance                                                             | Percentage of coating weight builds up                                             | Appearance (OK / Not OK) | Checked By Production Sign & Date | Verified By QA Sign & Date |
|------------------------------------------------------------------------|------------------------------------------------------------------------------------|--------------------------|-----------------------------------|----------------------------|
| White, Round, Biconvex, Film Coated tablets having scored on one side. | [Avg. wt. of coated tablet – Avg. wt. of core tablets] / Avg. wt. of core tablets. |                          |                                   |                            |

**8.9 Calculation of % Coating Weight Build up:**

| Avg. wt. of core 100 Tablets | Avg. wt. of 100 Coated tablet | % Of Coating Wt. Build up (Limit: 1.0 to 3.0 % w/w) | Checked By Production Sign & Date | Verified By QA Sign & Date |
|------------------------------|-------------------------------|-----------------------------------------------------|-----------------------------------|----------------------------|
|                              |                               |                                                     |                                   |                            |

**8.10 In-Process Specifications:**

| Sr. No. | Tests                          | Specifications                                                         | Limits                                                                 | Frequency             | Sample Qty. |
|---------|--------------------------------|------------------------------------------------------------------------|------------------------------------------------------------------------|-----------------------|-------------|
| 1.      | Description                    | White, Round, Biconvex, Film Coated tablets having scored on one side. | White, Round, Biconvex, Film Coated tablets having scored on one side. | At the end of coating | NA          |
| 2.      | Theoretical Weight of 20 tabs. | To comply                                                              | 2.660gm $\pm$ 2.0% (2.607gm to 2.713gm)                                | At the end of coating | 20 tabs     |
| 3.      | Uniformity of wt.              | To comply                                                              | 133.00mg $\pm$ 7.5% (123.025mg to 142.975 mg)                          | At the end of coating | 20 tabs     |
| 4.      | Thickness                      | To comply                                                              | 3.25 $\pm$ 0.2mm (3.05mm - 3.45mm)                                     | At the end of coating | 10 tabs     |
| 5.      | Disintegration time            | To comply                                                              | NMT 30 min                                                             | At the end of coating | 6 tabs      |



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date: 16/06/2020**

**Batch No.:**

**Batch Size:**



**8.11 In process Observation:**

| Parameter                      | Observation By |    |
|--------------------------------|----------------|----|
|                                | Production     | QA |
| Description                    |                |    |
| Theoretical Weight of 20 tabs. |                |    |
| Uniformity of wt               |                |    |
| Diameter (mm)                  |                |    |
| Thickness (mm)                 |                |    |
| DT (mins)                      |                |    |

**8.12 Coated Tablet Weight:**

| Container No.           | Clean status<br>OK/Not OK | Gross Wt.<br>(Kg) | Tare Wt.<br>(Kg) | Net Wt.<br>(Kg) | Remarks with<br>Sign & Date<br>(If any) |
|-------------------------|---------------------------|-------------------|------------------|-----------------|-----------------------------------------|
|                         |                           |                   |                  |                 |                                         |
|                         |                           |                   |                  |                 |                                         |
|                         |                           |                   |                  |                 |                                         |
|                         |                           |                   |                  |                 |                                         |
|                         |                           |                   |                  |                 |                                         |
|                         |                           |                   |                  |                 |                                         |
|                         |                           |                   |                  |                 |                                         |
|                         |                           |                   |                  |                 |                                         |
| Total No. of Container: |                           |                   | Total<br>Weight  |                 |                                         |

**Done By Production**  
**Sign & Date**

**Checked By Production**  
**Sign & Date**



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**8.13 Reconciliation of Coated Tablets:**

| S. No. | Calculation Points                      | Observed Values | Checked By<br>Production<br>Sign & Date | Verified By<br>QA<br>Sign & Date |
|--------|-----------------------------------------|-----------------|-----------------------------------------|----------------------------------|
| 1.     | Weight of Coated Tablets (Kg) (A)       |                 |                                         |                                  |
| 2.     | Theoretical Batch size (Kg) (B)         |                 |                                         |                                  |
| 3.     | In- Process Sample                      |                 |                                         |                                  |
| 4.     | QC Sample                               |                 |                                         |                                  |
| 5.     | Validation Sample                       |                 |                                         |                                  |
| 6.     | Total Sample (C)                        |                 |                                         |                                  |
| 7.     | Actual Yield $[(A + C / B) \times 100]$ |                 |                                         |                                  |
| 8.     | Non Recoverable (Kg)                    |                 |                                         |                                  |

**8.14 Sampling, Analysis and Release of Tablets:**

Send the intimation slip to QA Dept. to draw the sample of Tablets for analysis as per SOP No.: APD/QAD/028

**Sampling, Analysis Details**

| Intimated By<br>Production<br>Sign & Date | Intimation Received By<br>QA Sign & Date | Sampled By QA<br>QA Sign & Date | Quantity Sampled |
|-------------------------------------------|------------------------------------------|---------------------------------|------------------|
|                                           |                                          |                                 |                  |

The coated tablets are **Released/ Not Released** for Inspection/ Packing.

**QA Sign & Date:** \_\_\_\_\_

**Time:** \_\_\_\_\_

Attach the Report from QC: \_\_\_\_\_

**A.R. No:** \_\_\_\_\_

**8.15 Transfer above coated tablets to quarantine**

| Date of Transfer | Qty. Transferred<br>(in kg) | Transferred By<br>Sign & Date | Checked By QA<br>Sign & Date |
|------------------|-----------------------------|-------------------------------|------------------------------|
|                  |                             |                               |                              |

**8.16 Issuance of tablets for Inspection :**

| Date of Transfer | Date of Receiving | Qty. Received<br>in kg | Received By<br>Sign & Date | Verified By QA<br>Sign & Date |
|------------------|-------------------|------------------------|----------------------------|-------------------------------|
|                  |                   |                        |                            |                               |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:** 16/06/2020**Batch No.:****Batch Size:****8.17 Inspection Of Tablets****8.17.1 Line Clearance**

Area : ..... Line clearance checked at ..... on.....

Previous Product: ..... B. No.....

If the previous product same ensure that the powder of earlier batch is completely removed

| S. No.                         | Particulars                      | Remark (Yes/No) | Done By | Checked By Prod. | Verified By QA |
|--------------------------------|----------------------------------|-----------------|---------|------------------|----------------|
| 1.                             | Cleanliness of Area              |                 |         |                  |                |
| 2.                             | Cleanliness of Inspection belt   |                 |         |                  |                |
| 3.                             | Status Label of previous product |                 |         |                  |                |
| 4.                             | Container of previous product    |                 |         |                  |                |
| 5.                             | Tablet of previous product       |                 |         |                  |                |
| Temperature of area (NMT 25°C) |                                  |                 |         | °C               |                |
| RH of area (NMT 55%)           |                                  |                 |         | %                |                |

**Checked By Production Sign/Date..... Verified By QA Sign/Date .....****8.17.2 Tablet Verification:**

| Date | Time |    | Weight of tablets taken for Checking (kg) | Weight of good Tablets (kg) | Weight of non-recoverable residues (kg) | Weight of non-recoverable residues (kg) | Done by | Ckd. by |
|------|------|----|-------------------------------------------|-----------------------------|-----------------------------------------|-----------------------------------------|---------|---------|
|      | From | To |                                           |                             |                                         |                                         |         |         |
|      |      |    |                                           |                             |                                         |                                         |         |         |
|      |      |    |                                           |                             |                                         |                                         |         |         |
|      |      |    |                                           |                             |                                         |                                         |         |         |
|      |      |    |                                           |                             |                                         |                                         |         |         |

**8.17.3 Weight Verification:**

| Date         | Container No. | Gross weight (kg) | Tare weight (kg) | Net weight (kg) | Checked by |    |
|--------------|---------------|-------------------|------------------|-----------------|------------|----|
|              |               |                   |                  |                 | Prod.      | QA |
|              |               |                   |                  |                 |            |    |
|              |               |                   |                  |                 |            |    |
|              |               |                   |                  |                 |            |    |
|              |               |                   |                  |                 |            |    |
| Total weight |               |                   |                  |                 |            |    |



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**8.17.4 Finished Tablet Release:**

**Finished Tablet Release report attached here.**

**Report No. ....**

**Checked By Production Sign/Date..... Verified By QA Sign/Date .....**

**8.17.5 Yield Reconciliation of Inspected tablet :**

|                                     |   |                          |    |            |
|-------------------------------------|---|--------------------------|----|------------|
| Batch size (A)                      | : |                          | Kg | Tablets    |
| Tablets transferred for packing (B) | : |                          | Kg | Tablets    |
| Recoverable residues (C)            | : |                          | Kg | Tablets    |
| Non-recoverable residues (D)        | : |                          | Kg | Tablets    |
| Sample to QA analysis (E)           | : |                          |    | Tablets    |
| QC Sample (F)                       | : |                          |    | Tablets    |
| Total (G) = (B + C + E + F)         | : |                          | Kg | Tablets    |
| % Yield.                            | : | $\frac{G}{A} \times 100$ |    | -----<br>% |

**Limit: NLT 95.0%,**

**If yield is not within the limit, Justify:**


**Checked By Production Sign/Date..... Verified By QA Sign/Date .....**

**8.18 Destruction of Non Recoverable: Carry out Destruction as per SOP No.: APD/PRD/026.**

| Stage                | Non Recoverable Qty.(kg) | Qty. Destroyed (kg) | Destroyed By Sign& Date |
|----------------------|--------------------------|---------------------|-------------------------|
| Blending/Lubrication |                          |                     |                         |
| Compression          |                          |                     |                         |
| Coating              |                          |                     |                         |
| Inspection           |                          |                     |                         |

**Destroyed in presence of QA**

**Sign & Date: \_\_\_\_\_**

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:** 16/06/2020**Batch No.:****Batch Size:****9.0 ACCOUNTABILITY:**

| Sr. No. | Particulars                                      | In Kg | In No. |
|---------|--------------------------------------------------|-------|--------|
| 1.      | Theoretical batch quantity                       |       |        |
| 2.      | Total quantity transferred to packing department |       |        |
| 3.      | Blending loss                                    |       |        |
| 4.      | Compression loss                                 |       |        |
| 5.      | Coating loss                                     |       |        |
| 6.      | In process Sample                                |       |        |
| 7.      | Analysis Sample                                  |       |        |
| 8.      | Validation sample                                |       |        |

Accountability =  $\frac{(2+3 + 4+ 5 + 6+7+8)}{1} \times 100$

( ) X 100

= -----% (To be ascertained)

**Officer/ Executive Prod.**  
Sign and Date

**Head Production**  
Sign and Date

**Verified By QA**  
Sign and Date

**10.0 SIGNATURE LOG:**

| Sr. No. | Name of the Employee | Employee Code | Department | Specimen Signature |
|---------|----------------------|---------------|------------|--------------------|
| 1       |                      |               |            |                    |
| 2       |                      |               |            |                    |
| 3       |                      |               |            |                    |
| 4       |                      |               |            |                    |
| 5       |                      |               |            |                    |
| 6       |                      |               |            |                    |
| 7       |                      |               |            |                    |
| 8       |                      |               |            |                    |
| 9       |                      |               |            |                    |
| 10      |                      |               |            |                    |

**11.0 INCIDENT/DEVIATION RECEIVED (If Any):****Incident No.:****Deviation No:****Status:****Checked By QA:****Sign & Date:**



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name:** ALLYLESTRENOL TABLETS

**BMR No.:** APD/BMR/HT/ALL/001

**Revision No.:** 05

**BMR Supersedes No.:** 04

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



**15.0 BREAK DOWN DETAILS:**

| Date | From | To | Nature of Break Down | Area | Stage | Break Down Time | Rectification Time |
|------|------|----|----------------------|------|-------|-----------------|--------------------|
|      |      |    |                      |      |       |                 |                    |
|      |      |    |                      |      |       |                 |                    |
|      |      |    |                      |      |       |                 |                    |

**16.0 BACTH MANUFACTURING HISTORY SHEET:**

| S. No. | Stage         | Start Date | End Date | Checked By Production Sign & Date | Verified By QA Sign & Date |
|--------|---------------|------------|----------|-----------------------------------|----------------------------|
| 1.     | Dispensing    |            |          |                                   |                            |
| 2.     | Manufacturing |            |          |                                   |                            |
| 3.     | Compression   |            |          |                                   |                            |
| 4.     | Inspection    |            |          |                                   |                            |

**17. CHECK LIST FOR BATCH MANUFACTURING RECORDS & RELEASE FOR PACKING:**

| Sr. No. | Check Points                                                                                                                   | Status   |
|---------|--------------------------------------------------------------------------------------------------------------------------------|----------|
| 1.      | All the pages are available and comply with index.                                                                             | Yes / No |
| 2.      | Batch number available in all pages.                                                                                           | Yes / No |
| 3.      | Manufacturing and expiry date is correctly quoted.                                                                             | Yes / No |
| 4.      | Dispensing carried out on calibrated balance.                                                                                  | Yes / No |
| 5.      | Line Clearance is taken prior to all dispensing, Blending Area, Compression area activity.                                     | Yes / No |
| 6.      | Dispensing is carried out as per Bill of materials.                                                                            | Yes / No |
| 7.      | The entire dispensed label is properly affixed.                                                                                | Yes / No |
| 8.      | Environmental condition complies during all the manufacturing steps.                                                           | Yes / No |
| 9.      | No overwriting are observed all wrong entries are having been strike out and signed.                                           | Yes / No |
| 10.     | Actual equipments are used as specified.                                                                                       | Yes / No |
| 11.     | Batch dispensing is carried out after QA approval.                                                                             | Yes / No |
| 12.     | Batch manufacturing is carried out after QA approval.                                                                          | Yes / No |
| 13.     | Batch compression is carried out after QA approval.                                                                            | Yes / No |
| 14.     | QC approval is available on test request form                                                                                  | Yes / No |
| 15.     | Yield and reconciliation of all stage is within acceptance limit.                                                              | Yes / No |
| 16.     | All the in process check result comply within the acceptance limit                                                             | Yes / No |
| 17.     | All the blank space is filled correctly without pending entries.                                                               | Yes / No |
| 18.     | Batch reconciliation is completed and complies within the limit (In case it is outside limit justification is available)       | Yes / No |
| 19.     | If there is any deviation during the in process, it is recorded in process History sheet and deviation is raised and approved. | Yes / No |
| 20.     | Status Labels of Machinery                                                                                                     | Yes / No |
| 21.     | TR of wash water, blend and TR of compressed tablets attached.                                                                 | Yes / No |
| 22.     | RM excess issue order (if applicable)                                                                                          | Yes / No |

**ACCENT PHARMACEUTICALS & DIAGNOSTICS**

FOREST ROAD, SOLAN, H.P. (INDIA)

**BATCH MANUFACTURING RECORD****Product Name: ALLYLESTRENOL TABLETS****BMR No.: APD/BMR/HT/ALL/001****Revision No.: 05****BMR Supersedes No.: 04****Effective Date:** 16/06/2020**Batch No.:****Batch Size:**

|     |                                 |          |
|-----|---------------------------------|----------|
| 23. | Deviation and its Justification | Yes / No |
| 24. | Signature of Authorized Persons | Yes / No |
| 25. | Legibility of contents          | Yes / No |

**BMR Checked By:**  
**Production (Sign & Date)****BMR Reviewed By:**  
**QA (Sign & Date)****18. QC Data review Checklist:**

| Sr. | Check points                                                                   | Status   |
|-----|--------------------------------------------------------------------------------|----------|
| 1.  | Analytical report complete and approved                                        | Yes / No |
| 2.  | Records are readable, unequivocal & correct                                    | Yes / No |
| 3.  | Test according to current authorized specification and standard test procedure | Yes / No |
| 4.  | Test results are according to approved specification                           | Yes / No |
| 5.  | All raw data along with signed copy of audit Trail                             | Yes / No |

**Reviewed By QA/QC**  
**Sign & Date****19. ABBREVIATION:**

|                    |                                  |
|--------------------|----------------------------------|
| %                  | : Percentage                     |
| #                  | : Mesh                           |
| °C                 | : Degree Centigrade              |
| &                  | : And                            |
| Approx.            | : Approximately                  |
| A.R.               | : Analytical Report              |
| Avg.               | : Average                        |
| B. No.             | : Batch Number                   |
| Bal.               | : Balance                        |
| BMR                | : Batch Manufacturing Record     |
| Ckd.               | : Checked                        |
| Dt.                | : Date                           |
| Exp.               | : Expiry                         |
| Gm/g               | : Gram                           |
| GMP                | : Good Manufacturing Practice    |
| IP                 | : Indian Pharmacopoeia           |
| kg                 | : kilogram                       |
| L                  | : Liter                          |
| L.A.               | : Labeled Amount                 |
| L.C.               | : Label Claim                    |
| LOD                | : Loss on Drying                 |
| Mfg.               | : Manufacturing                  |
| mg                 | : Milligram                      |
| min                | : Minutes                        |
| Kg/cm <sup>2</sup> | : Kilogram per centimeter square |
| mm                 | : Millimeter                     |
| NA                 | : Not Applicable                 |
| NLT                | : Not Less Than                  |
| NMT                | : Not More Than                  |



**ACCENT PHARMACEUTICALS & DIAGNOSTICS**  
**FOREST ROAD, SOLAN, H.P. (INDIA)**  
**BATCH MANUFACTURING RECORD**

**Product Name: ALLYLESTRENOL TABLETS**

**BMR No.: APD/BMR/HT/ALL/001**

**Revision No.: 05**

**BMR Supersedes No.: 04**

**Effective Date:** 16/06/2020

**Batch No.:**

**Batch Size:**



No. / Nos. : Number / Numbers  
Prod. : Production  
QA : Quality Assurance  
QC : Quality Control  
Qty. : Quantity  
q.s. : Quantity Sufficient  
Ref. : Reference  
RH : Relative Humidity  
RMG : Rapid Mixer Granulator  
RPM : Rotation Per Minute  
R.O. : Role on  
Sr. No. : Serial Number  
SOP : Standard Operating Procedure  
S.S : Stainless Steel  
TAB : Tablet  
Temp. : Temperature  
Wt. : Weight  
w/w : Weight by Weight  
± : Pluses and Minus  
DT : Disintegration Test  
TR : Test Report

**20. DOCUMENT REVISION HISTORY:**

| Revision No. | Effective Date | CCF No.        | Details of Revision                                                 | Reason For Revision                                         |
|--------------|----------------|----------------|---------------------------------------------------------------------|-------------------------------------------------------------|
| 00           | 01 APR 2015    | NA             | New BMR                                                             | New BMR                                                     |
| 01           | 08 JULY 2016   | APD/CC/1506    | Separate SOP for Preparation of BMR and BPR is Incorporated.        | Separate SOP for Preparation of BMR and BPR is Incorporated |
| 02           | 08 JULY 2018   | NA             | Periodic Review                                                     | Periodic Review                                             |
| 03           | 24 MAR 2019    | APD/CC/1919    | Re formulation of BMR as per SOP No APD/QAD/023                     | System up-gradation                                         |
| 04           | 04 MAY 2020    | APD/PRD/20-014 | Removal of Doshion P 544 D and Additional of Sodium Lauryl Sulphate | Change in formulation                                       |
| 05           | 16 JUNE 2020   | APD/PRD/20-023 | Change in formulation                                               | Change in formulation                                       |